


Ha
Higgins
Analytical



2006

HPLC Columns and Cartridges

Higgins Analytical HPLC column product part numbers are based on a 2-4-4 character pattern. For example:

TS-1521-C185

is the part number for a conventional 15cm x 2.1mm stainless steel column packed with Targa C18 5µm sorbent.

The two character prefix indicates *sorbent and hardware* type.

Prefix Character 1 (Sorbent)

X-XXXX-XXXX

- C = CLYPEUS
- H = HAISIL 100 or HAISIL 300
- K = Kromasil
- L = HAISIL HL
- P = PHALANX
- R = PROTEO300
- T = TARGA
- V = HEAVY

Prefix Character 2 (Hardware)

X-XXXX-XXXX

- C = Capellini capillary column
- S = Stainless steel HPLC Column
- P = PEEK HPLC column
- F = Stainless steel cartridge column
- K = PEEK cartridge column
- R = Sprite Direct-Connect column

The four center characters indicate the column or cartridge *length (cm) and internal diameter (mm)*.

Center Characters 1 and 2 (Length)

XX-XX-XXXX

- 03 = 3cm
- 04 = 4cm
- 05 = 5cm
- 07 = 7.5cm
- 10 = 10cm
- 12 = 12.5cm
- 15 = 15cm
- 20 = 20cm
- 22 = 22cm
- 25 = 25cm
- 30 = 30cm
- 50 = 50cm

The four center characters indicate the column or cartridge *length (cm) and internal diameter (mm)*.

Center Characters 3 and 4 (I.D.)

XX-XX-XX-XXXX

- 07 = 0.075mm
- 15 = 0.15mm
- M3 = 0.3mm
- M5 = 0.5mm
- M8 = 0.8mm
- 01 = 1.0mm
- 21 = 2.1mm
- 30 = 3.0mm
- 32 = 3.2mm
- 40 = 4.0mm
- 46 = 4.6mm
- 10 = 10mm
- 20 = 20mm

The four character suffix indicates *sorbent chemistry and particle size*.

Suffix Characters 1 to 3 (Phase)

XX-XXXX-XXXXX

- C04 = Butyl or C4
- C08 = Octyl or C8
- C18 = ODS or C18
- CNP = Cyano or Nitrile
- M08 = HAISIL 100 C8
- M18 = HAISIL 100 C18
- NH2 = Amino
- PHN = Phenyl
- SIL = Silica Gel
- W04 = 300Å C4
- W08 = 300Å C8
- W18 = 300Å C18

Suffix Character 4 (Particle Size)

XX-XXXX-XXXXX

- 1 = 10µm
- 2 = 10 - 15µm
- 3 = 3 or 3.5µm
- 5 = 4.5 or 5µm
- 7 = 15 - 20µm

Notes:

- (1) Standard configurations are highlighted.
- (2) Polymeric phases have a 3-4-4 pattern.

PLx-xxxx-xxxx

Cartridge Holders

Guards

Two types of holders for 2cm long PEEK and s.s. cartridges are offered.

- HK-GUARD-FF Female/Female
- HK-GUARD-FM Direct Connect

Analytical Cartridge Holders

Three internal diameters (2.1, 3.2, and 4.6mm) and three lengths are offered.

- HK-HOLD-40 for 4cm long cartridge
- HK-HOLD-50 for 5cm long cartridge
- HK-HOLD-100 for 10cm long cartridge

HAISpeed 10 Cartridge Holders

These 10mm i.d. cartridges are offered in three lengths.

- 10-HOLD-30 for 3cm long cartridge
- 10-HOLD-50 for 5cm long cartridge
- 10-HOLD-100 for 10cm long cartridge

HAISpeed 20 Cartridge Holder

20mm i.d. prep guard cartridges require a 25mm long cartridge holder.

- 20-HOLD-25 for 20 x 25mm cartridge



CLYPEUS A “Small Shield” for a Big Job

Phases	C8, C18, Cyano, Phenyl
Particle Sizes	4.5 and 10µm
Pore Size	120Å
Pore Volume	0.8mL/gm
Surface Area	330m ² /gm
Carbon%(w/w)	C18 = 18%
Phase type	Extremely well end-capped monofunctional phases
Silica Class	Type B

Applications

The high purity base silica and thoroughly end capped bonded phases insure that CLYPEUS columns and cartridges will perform well in demanding situations where many other column brands fail. Even polar compounds elute as sharp and symmetrical peaks permitting surprisingly short columns to be used for high speed analysis. In addition to C8 and C18 phases, CLYPEUS columns are available with phenyl or cyano functionality.



HAISIL 100 The Workhorse

Phases	Silica, C8, and C18
Particle Sizes	3µm and 5µm
Pore Size	100Å
Pore Volume	0.5mL/gm
Surface Area	190m ² /gm
Carbon%(w/w)	C8 = 7%(w/w), C18 = 12%(w/w)
Phase type	Monofunctional & fully endcapped
Silica Class	Type A

Applications

HAISIL 100 columns and cartridges are very good general purpose columns that are particularly well suited for small molecules and peptide applications. The HAISIL 100 column and cartridge efficiencies are unrivaled yet they represent one of the most economical product lines in the market today. In formats ranging from 0.17 to 20mm ID and from 20 to 250mm in length,



HAISIL 300 For Biomolecules

Phases	C4, C8, and C18
Particle Sizes	5 and 10µm
Pore Size	300Å
Pore Volume	0.7mL/gm
Surface Area	80m ² /gm
%Carbon (w/w)	C4 = 3%, C8 = 4%, C18 = 6%
Phase type	Monofunctional and fully endcapped
Silica Class	Type A

Applications

HAISIL 300 columns and cartridges are ideal for biomolecule separations such as peptides and small proteins. Optimized phase bonding and column packing techniques assure unsurpassed performance and column-to-column reproducibility. Short cartridge columns are economical and are ideally suited for LC/MS applications. The column and cartridge offering ranges from 150µm to 20mm I.D. which means that bioseparations can be optimized by selecting a format that best fits your analytical requirements.



HAISIL HL High Load

Phases	C18
Particle Sizes	3 and 5µm
Pore Size	100Å
Pore Volume	0.8mL/gm
Surface Area	300m ² /gm
%Carbon (w/w)	C18 = 18%
Phase type	Monofunctional & fully endcapped
Silica Class	Type B

Applications

HAISIL HL “high load” columns and cartridges are intended for applications requiring maximum loading and selectivity. The high surface area permits high bonded phase loading yet the monofunctional bonded phase and the absence of micro-pores eliminates the problems that plague many of the so-called “high capacity” columns currently on the market. The Type B silica backbone is pure, robust, and stable under harsh conditions. Users report satisfactory lifetime when used for prolong periods at pH as high as 9.



HEAVY Very Hydrophobic Loading

Phases	C18
Particle Sizes	4.5µm
Pore Size	100Å
Pore Volume	1.1mL/gm
Surface Area	450m ² /gm
%Carbon (w/w)	C18 = 22%
Phase type	Monofunctional & fully endcapped
Silica Class	Type B

New!

Applications

Keeping in mind that relative retention (α) is a much more powerful term than plates (N) is the expression for resolution in chromatography, choosing a column with different selectivity is almost always more helpful than just going to a smaller particle size when trying to improve a difficult separation. Our new HEAVY columns provide the highest retention of any column we manufacture.

PROTEO300 Robust Bioseparations




Phases	C4 and C18
Particle Sizes	4.5 and 10µm
Pore Size	300Å
Pore Volume	0.9mL/gm
Surface Area	100m ² /gm
%Carbon (w/w)	C4 = 3%, C18= 8%
Phase type	Monofunctional & fully endcapped
Silica Class	Type B

New!

Applications

Low pH operating conditions and high pH column washing steps are typical bioprocessing operating conditions that can severely limit the lifetime and performance of widepore columns typically used for protein and peptide analysis. Higgins Analytical's new PROTEO300 has high capacity for enhanced resolution and is robust enough to withstand prolonged use at extreme pH (1.5 - 10).

PHALANX Ultra High Density End Capping



Phases	C18
Particle Sizes	4.5µm
Pore Size	100Å
Pore Volume	1.1mL/gm
Surface Area	430m ² /gm
%Carbon (w/w)	C18 = 18%
Phase type	Monofunctional & fully endcapped
Silica Class	Type B

Applications

Phalanx is uniquely suitable for the efficient analysis of both acidic and basic analytes. In addition, its high capacity, good stability at high pH, and fast equilibration at extreme ranges of hydrophobicity make Phalanx a powerful option for challenging separation problems.

TARGA A Uniquely Shielded Surface



Phases	C8, C18
Particle Sizes	3, 4.5 and 10µm
Pore Size	120Å
Pore Volume	0.8mL/gm
Surface Area	330m ² /gm
Carbon%(w/w)	C18 = 18%
Phase type	Monofunctional with a unique polar and bulky end-capping
Silica Class	Type B

Applications

TARGA columns and cartridges are ideal for applications requiring an absence of secondary solute-silanol interaction. Basic and acidic drugs and peptides can be chromatographed efficiently with little or no buffer. Targa columns are particularly well suited for LC/MS applications where acetic acid or formic acid can be substituted for TFA resulting in significant increase in signal. While TFA is not required, it does not preclude its use on TARGA columns. TARGA exhibits extremely fast equilibration kinetics from high to low organic phase compositions. It works well in 100% aqueous conditions for hydrophilic applications also.

KROMASIL 100



Phases	Silica, C4, C8, and C18
Particle Sizes	3.5 and 5µm
Pore Size	100Å
Pore Volume	0.9mL/gm
Surface Area	330m ² /gm
Carbon%(w/w)	C8 = 12%(w/w), C18 = 20%(w/w)
Phase type	Monofunctional & fully endcapped
Silica Class	Type B

Applications

The only commercially available sorbent offered in our product family, Kromasil has become well known for its unique selectivity, robust mechanical and chemical structure, and wide use in pharmaceutical process scale facilities. Higgins Analytical's ability to manufacture Kromasil HPLC columns that have unrivaled performance has made us an economically competitive and routine supplier of very high quality Kromasil columns to laboratories and OEM customers on five continents.

SPRITE

DIRECT CONNECT HPLC COLUMN



The Sprite family of narrowbore HPLC columns is a unique high performance product for LC-MS and other high-speed, small volume applications. All three sizes of 2.1mm I.D. Sprite HPLC columns are individually tested. A Certificate of Analysis accompanies each column. As pictured on the back page, the Sprite's certificate includes the original chromatogram and comprehensive testing and sorbent documentation.

The inert PEEK and titanium design has been proven to be structurally sound to over 5,000 PSI (340Bar).

PICCOLO

SHORT MICROBORE COLUMNS



The Piccolo family of microbore HPLC columns is a unique high performance product for LC-MS and other high-speed, small volume applications. These efficient columns are the only products of their type that are individually tested. A Certificate of Analysis accompanies each column.

In addition to analytical microbore applications Piccolo columns are ideal guards for protecting longer microbore HPLC columns as illustrated on the back page.

Piccolo column void volumes range from approximately 5µl to 300nL, making them perfect for a wide variety of trapping, desalting, and trace enrichment applications. These small columns can be conveniently used with micro injection and switching valves as pictured below.

HAIPEEK CARTRIDGES

PEEK/TITANIUM CARTRIDGE COLUMNS



The HAIPEEK family of narrowbore HPLC cartridge columns is a unique high performance product for LC-MS and other high-speed, small volume applications. The 40 and 50mm sizes of 2.1mm I.D. HAIPEEK HPLC Cartridges are individually tested. A Certificate of Analysis, which includes the original chromatogram and comprehensive testing and sorbent documentation, accompanies each column. The 20mm x 2.1mm cartridges excel as narrowbore guard columns, yet even though they are not individually tested, the majority are routinely used as very economical and efficient columns for high speed HPLC screening and LC-MS applications. The back page describes the holders required for Higgins Analytical's HAIPEEK Cartridge Columns.

CAPELLINI

CAPILLARY MICROBORE COLUMNS



Higgins Analytical Capellini column are robust and high performance capillary HPLC columns for LC-MS and other very small scale applications. The inert PEEK and titanium fittings complement the Higgins Analytical sorbents such as TARGA which was developed to enable LC-MS applications without TFA.

Our no nonsense testing procedure does not resort to sample and flow splitting techniques, inaccurate k' values, and miniscule chromatograms, as do some manufacturers, to illustrate column performance. The CofA's symmetrical peaks and the high resolution map speak for themselves.

HAI Speed 10 CARTRIDGES

10mm I.D. SemiPrep Cartridge Columns



HAI Speed 10 - A high performance 10mm I.D. Cartridge System.

High throughput requirements for combinatorial chemistry applications often call for 10mm I.D. columns. Higgins Analytical has developed a convenient cartridge system that offers an economical solution for demanding combinatorial and semipreparative needs. Even though a HAI Speed 10 cartridge costs a fraction of conventional 10mm HPLC column of the same length, they have the same consistently high performance of other Higgins Analytical columns.

HAI Speed 20 CARTRIDGES

20mm I.D. Prep Guard Cartridge Columns




HAI Speed 20 guard cartridges provide unique Protection for 20mm I.D. prep columns.

"Prep is a dirty business," and expensive columns can be quickly destroyed by crude samples if proper precautions are not taken. Higgins Analytical's HAI Speed 20 Prep Guard System provides effective column protection without compromising performance. In addition to an efficiently packed bed which retains strongly absorbing junk that would change the selectivity of the prep column, these unique guard cartridges have 0.5µm frits that remove particulate contaminants as well.

HAIPEEK & HAICART Guards

PEEK/TITANIUM and s.s. Guard Cartridges



The complete range of sorbents manufactured by Higgins Analytical is available in 2cm x 2.1 or 3.2mm I.D. guard cartridges. Not only does this provide a wide selection of products for protecting analytical and narrowbore HPLC columns, it also facilitates a wide variety of trace enrichment, column switching and ballistic gradient applications.

These 2cm long cartridges require a reusable holder and are compatible with Higgins Analytical's HAIPEEK Guard Holders as well as PE/Brownlee's NewGuard Holder.

Cartridge Holders

Holders for 2.1, 3.2, 4.6, 10, and 20mm I.D. Cartridges



Even though Higgins Analytical cartridge holders do not require any tools for assembly and disassembly, their leak-free operation assures that maximum performance will be achieved in your demanding, high-speed separations.

The Holder's light weight also makes it easy to handle and assures its compatibility in challenging plumbing situations.

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
10-HOLD-100	HAI Speed 10 100mm Cartridge Holder Kit	
10-HOLD-30	HAI Speed 10 30mm Cartridge Holder Kit	
10-HOLD-50	HAI Speed 10 50mm Cartridge Holder Kit	
20-HOLD-25	HAI Speed 20 20x25mm Cartridge Holder Kit	
C108	CLIQUEUS C8 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
C118	CLIQUEUS C18 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
C1CN	CLIQUEUS CN 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
C1PH	CLIQUEUS PH 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
CC-0507-C185	CLIQUEUS Capellini C18 5µm, 50 x 0.075mm column	
CC-0515-C185	CLIQUEUS Capellini C18 5µm, 50 x 0.150mm column	
CC-1007-C185	CLIQUEUS Capellini C18 5µm, 100 x 0.075mm column	
CC-1015-C185	CLIQUEUS Capellini C18 5µm, 100 x 0.150mm column	
CC-1507-C185	CLIQUEUS Capellini C18 5µm, 150 x 0.075mm column	
CC-1515-C185	CLIQUEUS Capellini C18 5µm, 150 x 0.150mm column	
CC-2507-C185	CLIQUEUS Capellini C18 5µm, 250 x 0.075mm column	
CC-2515-C185	CLIQUEUS Capellini C18 5µm, 250 x 0.150mm column	
CF-0220-C181	CLIQUEUS C18 10µm, 25cm x 20mm Cartridge s.s. 2/pkg	
CF-0310-C085	CLIQUEUS C8 5µm, 3cm x 10mm Cartridge s.s. 1/pkg	
CF-0310-C181	CLIQUEUS C18 10µm, 3cm x 10mm Cartridge s.s. 1/pkg	
CF-0310-C185	CLIQUEUS C18 5µm, 3cm x 10mm Cartridge s.s. 1/pkg	
CF-0310-CNP5	CLIQUEUS Cyano 5µm, 3cm x 10mm Cartridge s.s. 1/pkg	
CF-0310-PHN5	CLIQUEUS Phenyl 5µm, 3cm x 10mm Cartridge s.s. 1/pkg	
CF-0432-C085	CLIQUEUS C8 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0432-C185	CLIQUEUS C18 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0432-CNP5	CLIQUEUS Cyano 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0432-PHN5	CLIQUEUS Phenyl 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0446-C085	CLIQUEUS C8 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
CF-0446-C185	CLIQUEUS C18 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
CF-0446-CNP5	CLIQUEUS Cyano 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
CF-0446-PHN5	CLIQUEUS Phenyl 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
CF-0510-C085	CLIQUEUS C8 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
CF-0510-C181	CLIQUEUS C18 10µm, 5cm x 10mm Cartridge s.s. 1/pkg	
CF-0510-C185	CLIQUEUS C18 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
CF-0510-CNP5	CLIQUEUS Cyano 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
CF-0510-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
CF-0532-C085	CLIQUEUS C8 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0532-C185	CLIQUEUS C18 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0532-CNP5	CLIQUEUS Cyano 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0532-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
CF-0546-C085	CLIQUEUS C8 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
CF-0546-C185	CLIQUEUS C18 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
CF-0546-CNP5	CLIQUEUS Cyano 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
CF-0546-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
CF-1010-C085	CLIQUEUS C8 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
CF-1010-C181	CLIQUEUS C18 10µm, 10cm x 10mm Cartridge s.s. 1/pkg	
CF-1010-C185	CLIQUEUS C18 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
CF-1010-CNP5	CLIQUEUS Cyano 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
CF-1010-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
CF-1032-C085	CLIQUEUS C8 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
CF-1032-C185	CLIQUEUS C18 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
CF-1032-CNP5	CLIQUEUS Cyano 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
CF-1032-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
CF-1046-C085	CLIQUEUS C8 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
CF-1046-C185	CLIQUEUS C18 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
CF-1046-CNP5	CLIQUEUS Cyano 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
CF-1046-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
CK-0221-C085	CLIQUEUS C8 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
CK-0221-C185	CLIQUEUS C18 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
CK-0221-CNP5	CLIQUEUS Cyano 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
CK-0221-PHN5	CLIQUEUS Phenyl 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
CK-0421-C085	CLIQUEUS C8 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
CK-0421-C185	CLIQUEUS C18 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
CK-0421-CNP5	CLIQUEUS Cyano 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
CK-0421-PHN5	CLIQUEUS Phenyl 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
CK-0521-C085	CLIQUEUS C8 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
CK-0521-C185	CLIQUEUS C18 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
CK-0521-CNP5	CLIQUEUS Cyano 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
CK-0521-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
CP-0101-C185	CLIQUEUS Piccolo C18 5µm, 10 x 1.0mm column	
CP-M2M5-C185	CLIQUEUS Piccolo C18 5µm, 2.5 x 0.5mm column	
CP-M501-C185	CLIQUEUS Piccolo C18 5µm, 5 x 1.0mm column	
CR-0121-C185	CLIQUEUS Sprite C18 5µm, 10 x 2.1mm column	
CR-0221-C185	CLIQUEUS Sprite C18 5µm, 20 x 2.1mm column	
CR-0421-C185	CLIQUEUS Sprite C18 5µm, 40 x 2.1mm column	
CS-0301-C085	CLIQUEUS C8 5µm, 3cm x 1.0mm Microbore Column s.s.	
CS-0301-C185	CLIQUEUS C18 5µm, 3cm x 1.0mm Microbore Column s.s.	
CS-0301-CNP5	CLIQUEUS Cyano 5µm, 3cm x 1.0mm Microbore Column s.s.	
CS-0301-PHN5	CLIQUEUS Phenyl 5µm, 3cm x 1.0mm Microbore Column s.s.	
CS-0321-C085	CLIQUEUS C8 5µm, 3cm x 2.1mm Analytical Column s.s.	
CS-0321-C185	CLIQUEUS C18 5µm, 3cm x 2.1mm Analytical Column s.s.	
CS-0321-CNP5	CLIQUEUS Cyano 5µm, 3cm x 2.1mm Analytical Column s.s.	
CS-0321-PHN5	CLIQUEUS Phenyl 5µm, 3cm x 2.1mm Analytical Column s.s.	
CS-0330-C085	CLIQUEUS C8 5µm, 3cm x 3.0mm Analytical Column s.s.	
CS-0330-C185	CLIQUEUS C18 5µm, 3cm x 3.0mm Analytical Column s.s.	
CS-0330-CNP5	CLIQUEUS Cyano 5µm, 3cm x 3.0mm Analytical Column s.s.	
CS-0330-PHN5	CLIQUEUS Phenyl 5µm, 3cm x 3.0mm Analytical Column s.s.	
CS-0346-C085	CLIQUEUS C8 5µm, 3cm x 4.6mm Analytical Column s.s.	
CS-0346-C181	CLIQUEUS C18 10µm, 3cm x 4.6mm Analytical Column s.s.	
CS-0346-C185	CLIQUEUS C18 5µm, 3cm x 4.6mm Analytical Column s.s.	
CS-0346-CNP5	CLIQUEUS Cyano 5µm, 3cm x 4.6mm Analytical Column s.s.	
CS-0346-PHN5	CLIQUEUS Phenyl 5µm, 3cm x 4.6mm Analytical Column s.s.	
CS-03M3-C085	CLIQUEUS C8 5µm, 3cm x 0.3mm Microbore Column s.s.	
CS-03M3-C185	CLIQUEUS C18 5µm, 3cm x 0.3mm Microbore Column s.s.	
CS-03M3-CNP5	CLIQUEUS Cyano 5µm, 3cm x 0.3mm Microbore Column s.s.	
CS-03M3-PHN5	CLIQUEUS Phenyl 5µm, 3cm x 0.3mm Microbore Column s.s.	
CS-03M5-C085	CLIQUEUS C8 5µm, 3cm x 0.5mm Microbore Column s.s.	
CS-03M5-C185	CLIQUEUS C18 5µm, 3cm x 0.5mm Microbore Column s.s.	
CS-03M5-CNP5	CLIQUEUS Cyano 5µm, 3cm x 0.5mm Microbore Column s.s.	
CS-03M5-PHN5	CLIQUEUS Phenyl 5µm, 3cm x 0.5mm Microbore Column s.s.	
CS-0501-C085	CLIQUEUS C8 5µm, 5cm x 1.0mm Microbore Column s.s.	
CS-0501-C185	CLIQUEUS C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
CS-0501-CNP5	CLIQUEUS Cyano 5µm, 5cm x 1.0mm Microbore Column s.s.	
CS-0501-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 1.0mm Microbore Column s.s.	
CS-0510-C085	CLIQUEUS C8 5µm, 5cm x 10mm Semi/Prep Column s.s.	
CS-0510-C181	CLIQUEUS C18 10µm 5cm x 10mm Semi/Prep Column s.s.	
CS-0510-C185	CLIQUEUS C18 5µm, 5cm x 10mm Semi/Prep Column s.s.	
CS-0520-C085	CLIQUEUS C8 5µm, 5cm x 20mm Semi/Prep Column s.s.	
CS-0520-C181	CLIQUEUS C18 10µm, 5cm x 20mm Semi/Prep Column s.s.	
CS-0520-C185	CLIQUEUS C18 5µm, 5cm x 20mm Semi/Prep Column s.s.	
CS-0521-C085	CLIQUEUS C8 5µm, 5cm x 2.1mm Analytical Column s.s.	
CS-0521-C185	CLIQUEUS C18 5µm, 5cm x 2.1mm Analytical Column s.s.	
CS-0521-CNP5	CLIQUEUS Cyano 5µm, 5cm x 2.1mm Analytical Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
CS-0521-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 2.1mm Analytical Column s.s.	
CS-0530-C085	CLIQUEUS C8 5µm, 5cm x 3.0mm Analytical Column s.s.	
CS-0530-C185	CLIQUEUS C18 5µm, 5cm x 3.0mm Analytical Column s.s.	
CS-0530-CNP5	CLIQUEUS Cyano 5µm, 5cm x 3.0mm Analytical Column s.s.	
CS-0530-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 3.0mm Analytical Column s.s.	
CS-0546-C085	CLIQUEUS C8 5µm, 5cm x 4.6mm Analytical Column s.s.	
CS-0546-C181	CLIQUEUS C18 10µm 5cm x 4.6mm Analytical Column s.s.	
CS-0546-C185	CLIQUEUS C18 5µm, 5cm x 4.6mm Analytical Column s.s.	
CS-0546-CNP5	CLIQUEUS Cyano 5µm, 5cm x 4.6mm Analytical Column s.s.	
CS-0546-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 4.6mm Analytical Column s.s.	
CS-05M3-C085	CLIQUEUS C8 5µm, 5cm x 0.3mm Microbore Column s.s.	
CS-05M3-C185	CLIQUEUS C18 5µm, 5cm x 0.3mm Microbore Column s.s.	
CS-05M3-CNP5	CLIQUEUS Cyano 5µm, 5cm x 0.3mm Microbore Column s.s.	
CS-05M3-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 0.3mm Microbore Column s.s.	
CS-05M5-C085	CLIQUEUS C8 5µm, 5cm x 0.5mm Microbore Column s.s.	
CS-05M5-C185	CLIQUEUS C18 5µm, 5cm x 0.5mm Microbore Column s.s.	
CS-05M5-CNP5	CLIQUEUS Cyano 5µm, 5cm x 0.5mm Microbore Column s.s.	
CS-05M5-PHN5	CLIQUEUS Phenyl 5µm, 5cm x 0.5mm Microbore Column s.s.	
CS-1001-C085	CLIQUEUS C8 5µm, 10cm x 1.0mm Microbore Column s.s.	
CS-1001-C185	CLIQUEUS C18 5µm, 10cm x 1.0mm Microbore Column s.s.	
CS-1001-CNP5	CLIQUEUS Cyano 5µm, 10cm x 1.0mm Microbore Column s.s.	
CS-1001-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 1.0mm Microbore Column s.s.	
CS-1010-C085	CLIQUEUS C8 5µm, 10cm x 10mm Semi/Prep Column s.s.	
CS-1010-C181	CLIQUEUS C18 10µm, 10cm x 10mm Semi/Prep Column s.s.	
CS-1010-C185	CLIQUEUS C18 5µm, 10cm x 10mm Semi/Prep Column s.s.	
CS-1020-C085	CLIQUEUS C8 5µm, 10cm x 20mm Semi/Prep Column s.s.	
CS-1020-C181	CLIQUEUS C18 10µm, 10cm x 20mm Semi/Prep Column s.s.	
CS-1020-C185	CLIQUEUS C18 5µm, 10cm x 20mm Semi/Prep Column s.s.	
CS-1021-C085	CLIQUEUS C8 5µm, 10cm x 2.1mm Analytical Column s.s.	
CS-1021-C185	CLIQUEUS C18 5µm, 10cm x 2.1mm Analytical Column s.s.	
CS-1021-CNP5	CLIQUEUS Cyano 5µm, 10cm x 2.1mm Analytical Column s.s.	
CS-1021-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 2.1mm Analytical Column s.s.	
CS-1030-C085	CLIQUEUS C8 5µm, 10cm x 3.0mm Analytical Column s.s.	
CS-1030-C185	CLIQUEUS C18 5µm, 10cm x 3.0mm Analytical Column s.s.	
CS-1030-CNP5	CLIQUEUS Cyano 5µm, 10cm x 3.0mm Analytical Column s.s.	
CS-1030-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 3.0mm Analytical Column s.s.	
CS-1046-C085	CLIQUEUS C8 5µm, 10cm x 4.6mm Analytical Column s.s.	
CS-1046-C181	CLIQUEUS C18 10µm, 10cm x 4.6mm Analytical Column s.s.	
CS-1046-C185	CLIQUEUS C18 5µm, 10cm x 4.6mm Analytical Column s.s.	
CS-1046-CNP5	CLIQUEUS Cyano 5µm, 10cm x 4.6mm Analytical Column s.s.	
CS-1046-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 4.6mm Analytical Column s.s.	
CS-10M3-C085	CLIQUEUS C8 5µm, 10cm x 0.3mm Microbore Column s.s.	
CS-10M3-C185	CLIQUEUS C18 5µm, 10cm x 0.3mm Microbore Column s.s.	
CS-10M3-CNP5	CLIQUEUS Cyano 5µm, 10cm x 0.3mm Microbore Column s.s.	
CS-10M3-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 0.3mm Microbore Column s.s.	
CS-10M5-C085	CLIQUEUS C8 5µm, 10cm x 0.5mm Microbore Column s.s.	
CS-10M5-C185	CLIQUEUS C18 5µm, 10cm x 0.5mm Microbore Column s.s.	
CS-10M5-CNP5	CLIQUEUS Cyano 5µm, 10cm x 0.5mm Microbore Column s.s.	
CS-10M5-PHN5	CLIQUEUS Phenyl 5µm, 10cm x 0.5mm Microbore Column s.s.	
CS-1501-C085	CLIQUEUS C8 5µm, 15cm x 1.0mm Microbore Column s.s.	
CS-1501-C185	CLIQUEUS C18 5µm, 15cm x 1.0mm Microbore Column s.s.	
CS-1501-CNP5	CLIQUEUS Cyano 5µm, 15cm x 1.0mm Microbore Column s.s.	
CS-1501-PHN5	CLIQUEUS Phenyl 5µm, 15cm x 1.0mm Microbore Column s.s.	
CS-1510-C085	CLIQUEUS C8 5µm, 15cm x 10mm Semi/Prep Column s.s.	
CS-1510-C181	CLIQUEUS C18 10µm, 15cm x 10mm Semi/Prep Column s.s.	
CS-1510-C185	CLIQUEUS C18 5µm, 15cm x 10mm Semi/Prep Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
CS-1520-C085	CLIQUEUS C8 5µm, 15cm x 20mm Semi/Prep Column s.s.	
CS-1520-C181	CLIQUEUS C18 10µm, 15cm x 20mm Semi/Prep Column s.s.	
CS-1520-C185	CLIQUEUS C18 5µm, 15cm x 20mm Semi/Prep Column s.s.	
CS-1521-C085	CLIQUEUS C8 5µm, 15cm x 2.1mm Analytical Column s.s.	
CS-1521-C185	CLIQUEUS C18 5µm, 15cm x 2.1mm Analytical Column s.s.	
CS-1521-CNP5	CLIQUEUS Cyano 5µm, 15cm x 2.1mm Analytical Column s.s.	
CS-1521-PHN5	CLIQUEUS Phenyl 5µm, 15cm x 2.1mm Analytical Column s.s.	
CS-1530-C085	CLIQUEUS C8 5µm, 15cm x 3.0mm Analytical Column s.s.	
CS-1530-C185	CLIQUEUS C18 5µm, 15cm x 3.0mm Analytical Column s.s.	
CS-1530-CNP5	CLIQUEUS Cyano 5µm, 15cm x 3.0mm Analytical Column s.s.	
CS-1530-PHN5	CLIQUEUS Phenyl 5µm, 15cm x 3.0mm Analytical Column s.s.	
CS-1546-C085	CLIQUEUS C8 5µm, 15cm x 4.6mm Analytical Column s.s.	
CS-1546-C181	CLIQUEUS C18 10µm, 15cm x 4.6mm Analytical Column s.s.	
CS-1546-C185	CLIQUEUS C18 5µm, 15cm x 4.6mm Analytical Column s.s.	
CS-1546-CNP5	CLIQUEUS Cyano 5µm, 15cm x 4.6mm Analytical Column s.s.	
CS-1546-PHN5	CLIQUEUS Phenyl 5µm, 15cm x 4.6mm Analytical Column s.s.	
CS-15M5-C085	CLIQUEUS C8 5µm, 15cm x 0.5mm Microbore Column s.s.	
CS-15M5-C185	CLIQUEUS C18 5µm, 15cm x 0.5mm Microbore Column s.s.	
CS-15M5-CNP5	CLIQUEUS Cyano 5µm, 15cm x 0.5mm Microbore Column s.s.	
CS-15M5-PHN5	CLIQUEUS Phenyl 5µm, 15cm x 0.5mm Microbore Column s.s.	
CS-2501-C085	CLIQUEUS C8 5µm, 25cm x 1.0mm Microbore Column s.s.	
CS-2501-C185	CLIQUEUS C18 5µm, 25cm x 1.0mm Microbore Column s.s.	
CS-2501-CNP5	CLIQUEUS Cyano 5µm, 25cm x 1.0mm Microbore Column s.s.	
CS-2501-PHN5	CLIQUEUS Phenyl 5µm, 25cm x 1.0mm Microbore Column s.s.	
CS-2510-C085	CLIQUEUS C8 5µm, 25cm x 10mm Semi/Prep Column s.s.	
CS-2510-C181	CLIQUEUS C18 10µm, 25cm x 10mm Semi/Prep Column s.s.	
CS-2510-C185	CLIQUEUS C18 5µm, 25cm x 10mm Semi/Prep Column s.s.	
CS-2520-C085	CLIQUEUS C8 5µm, 25cm x 20mm Semi/Prep Column s.s.	
CS-2520-C181	CLIQUEUS C18 10µm, 25cm x 20mm Semi/Prep Column s.s.	
CS-2520-C185	CLIQUEUS C18 5µm, 25cm x 20mm Semi/Prep Column s.s.	
CS-2521-C085	CLIQUEUS C8 5µm, 25cm x 2.1mm Analytical Column s.s.	
CS-2521-C185	CLIQUEUS C18 5µm, 25cm x 2.1mm Analytical Column s.s.	
CS-2521-CNP5	CLIQUEUS Cyano 5µm, 25cm x 2.1mm Analytical Column s.s.	
CS-2521-PHN5	CLIQUEUS Phenyl 5µm, 25cm x 2.1mm Analytical Column s.s.	
CS-2530-C085	CLIQUEUS C8 5µm, 25cm x 3.0mm Analytical Column s.s.	
CS-2530-C185	CLIQUEUS C18 5µm, 25cm x 3.0mm Analytical Column s.s.	
CS-2530-CNP5	CLIQUEUS Cyano 5µm, 25cm x 3.0mm Analytical Column s.s.	
CS-2530-PHN5	CLIQUEUS Phenyl 5µm, 25cm x 3.0mm Analytical Column s.s.	
CS-2546-C085	CLIQUEUS C8 5µm, 25cm x 4.6mm Analytical Column s.s.	
CS-2546-C181	CLIQUEUS C18 10µm, 25cm x 4.6mm Analytical Column s.s.	
CS-2546-C185	CLIQUEUS C18 5µm, 25cm x 4.6mm Analytical Column s.s.	
CS-2546-CNP5	CLIQUEUS Cyano 5µm, 25cm x 4.6mm Analytical Column s.s.	
CS-2546-PHN5	CLIQUEUS Phenyl 5µm, 25cm x 4.6mm Analytical Column s.s.	
FF-FKIT	HAIFILTER FF (Female/Female) Starter Kit	
GD-AION	Anion Exchange, 0.2µm porosity 10/pkg	
GD-C18A	ODS (Optional Teflon, 0.2µm porosity 10/pkg	
GD-C18B	ODS (Sandard PVC, 0.8mm porosity)	
GD-CION	Cation Exchange, 0.2µm porosity 10/pkg	
GD-GKIT	HAIGUARD Holder (female/male) with GD-C18B discs (pkg of 10)	
GD-HOLD	HAIGUARD Holder (female/male) without discs	
GD-HOLW	GD-WKIT for Waters HPLC column (without Guard Disc	
GD-SDVB	Styrene divinylbenzene, 0.2µm porosity 10/pkg	
GD-WKIT	HAIGUARD Holder (female/male for Waters Columns) with GD-C18B discs (pkg of 10)	
GF-GKIT	HAIGUARD Holder (female/female) with GD-C18B discs (pkg of 10)	
GF-HOLD	HAIGUARD Holder (female/female) without discs	
HI08	HAISIL 100 C8 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
H118	HAISIL 100 C18 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
H151	HAISIL 100 Silica 5µm, 2cm 3.2mm Cartridge s.s. 5/pkg	
H304	HAISIL 300 C4 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
H308	HAISIL 300 C8 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
H318	HAISIL 300 C18 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
HF-0220-M181	HAISIL 100 C18 10µm, 2.5cm x 20mm Cartridge 2/pkg	
HF-0220-SIL1	HAISIL 100 Silica 10µm, 2.5cm x 20mm Cartridge 2/pkg	
HF-0220-W181	HAISIL 300 C18 10µm, 2.5cm x 20mm Cartridge s.s. 2/pkg	
HF-0220-W182	HAISIL 300 C18 12µm, 2.5 x 20mm Cartridge s.s. 2/pkg	
HF-0310-M085	HAISIL 100 C8 5µm, 3cm x 10mm s.s. Cartridge 1/pkg	
HF-0310-M185	HAISIL 100 C18 5µm, 3cm x 10mm s.s. Cartridge 1/pkg	
HF-0310-SIL5	HAISIL 100 Silica 5µm, 3cm x 10mm s.s. Cartridge 1/pkg	
HF-0310-W045	HAISIL 300 C4 5µm, 3cm x 10mm s.s. Cartridge 1/pkg	
HF-0310-W085	HAISIL 300 C8 5µm, 3cm x 10mm s.s. Cartridge 1/pkg	
HF-0310-W185	HAISIL 300 C18 5µm, 3cm x 10mm s.s. Cartridge 1/pkg	
HF-0432-M085	HAISIL 100 C8 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0432-M185	HAISIL 100 C18 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0432-SIL5	HAISIL 100 Silica 5µm, 4cm x 3.2mm s.s. Cartridge 2/pkg	
HF-0432-W045	HAISIL 300 C4 5µm, 4cm x 3.2mm s.s. Cartridge 2/pkg	
HF-0432-W085	HAISIL 300 C8 5µm, 4cm x 3.2mm s.s. Cartridge 2/pkg	
HF-0432-W185	HAISIL 300 C18 5µm, 4cm x 3.2mm s.s. Cartridge 2/pkg	
HF-0446-M085	HAISIL 100 C8 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0446-M185	HAISIL 100 C18 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0446-SIL5	HAISIL 100 Silica 5µm, 4cm x 4.6mm s.s. Cartridge 2/pkg	
HF-0446-W045	HAISIL 300 C4 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0446-W085	HAISIL 300 C8 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0446-W185	HAISIL 300 C18 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0510-M085	HAISIL 100 C8 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
HF-0510-M185	HAISIL 100 C18 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
HF-0510-SIL5	HAISIL 100 Silica 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
HF-0510-W045	HAISIL 300 C4 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
HF-0510-W085	HAISIL 300 C8 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
HF-0510-W185	HAISIL 300 C18 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
HF-0532-M085	HAISIL 100 C8 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0532-M185	HAISIL 100 C18 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0532-SIL5	HAISIL 100 Silica 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0532-W045	HAISIL 300 C4 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0532-W085	HAISIL 300 C8 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0532-W185	HAISIL 300 C18 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
HF-0546-M085	HAISIL 100 C8 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0546-M185	HAISIL 100 C18 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0546-SIL5	HAISIL 100 Silica 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0546-W045	HAISIL 300 C4 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0546-W085	HAISIL 300 C8 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
HF-0546-W185	HAISIL 300 C18 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
HF-1010-M085	HAISIL 100 C18 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
HF-1010-M185	HAISIL 100 C18 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
HF-1010-SIL5	HAISIL 100 Silica 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
HF-1010-W045	HAISIL 300 C4 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
HF-1010-W085	HAISIL 300 C8 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
HF-1010-W185	HAISIL 300 C18 5µm, 10cm x 10mm Cartridge s.s. 1/pkg	
HF-1032-M085	HAISIL 100 C8 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
HF-1032-M185	HAISIL 100 C18 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
HF-1032-SIL5	HAISIL 100 Silica 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
HF-1032-W045	HAISIL 300 C4 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
HF-1032-W085	HAISIL 300 C8 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
HF-1032-W185	HAISIL 300 C18 5µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
HF-1046-M085	HAISIL 100 C8 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
HF-1046-M185	HAISIL 100 C18 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
HF-1046-SIL5	HAISIL 100 Silica 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
HF-1046-W045	HAISIL 300 C4 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
HF-1046-W085	HAISIL 300 C8 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
HF-1046-W185	HAISIL 300 C18 5µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
HF-FKIT	HAIFILTER FM (Standard)	
HF-PFILT	PEEK HAIFILTER replacements filters (pkg of 10)	
HF-RKIT	HAIFILTER FM (for Rheodyne injection valves -	
HF-SFILT	Stainless steel HAIFILTER replacement filters (pkg of 10)	
HF-TFILT	Titanium HAIFILTER replacements filters (pkg of 10)	
HF-WKIT	HAIFILTER FM (for Waters HPLC columns)	
HK-0221-M085	HAISIL 100 C8 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
HK-0221-M185	HAISIL 100 C18 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
HK-0221-SIL5	HAISIL 100 Silica 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
HK-0221-W045	HAISIL 300 C4 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
HK-0221-W085	HAISIL 300 C8 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
HK-0221-W185	HAISIL 300 C18 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
HK-0421-M085	HAISIL 100 C8 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0421-M185	HAISIL 100 C18 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0421-SIL5	HAISIL 100 Silica 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0421-W045	HAISIL 300 C4 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0421-W085	HAISIL 300 C8 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0421-W185	HAISIL 300 C18 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0521-M085	HAISIL 100 C8 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0521-M185	HAISIL 100 C18 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0521-SIL5	HAISIL 100 Silica 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0521-W045	HAISIL 300 C4 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0521-W085	HAISIL 300 C8 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
HK-0521-W185	HAISIL 300 C18 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
HK-GUARD-FF	HAIPEEK Guard Holder Kit (female/female)	
HK-GUARD-FM	HAIPEEK Guard Holder Kit (female/male)	
HK-HOLD-100	HAIPEEK 100mm Cartridge Holder Kit	
HK-HOLD-40	HAIPEEK 40mm Cartridge Holder Kit	
HK-HOLD-50	HAIPEEK 50mm Cartridge Holder Kit	
HP-0521-W045	HAISIL 300 C4 5µm, 5cm x 2.1mm Analytical Column PEEK	
HP-0521-W085	HAISIL 300 C8 5µm, 5cm x 2.1mm Analytical Column PEEK	
HP-0521-W185	HAISIL 300 C18 5µm, 5cm x 2.1mm Analytical Column PEEK	
HP-0530-W045	HAISIL 300 C4 5µm, 5cm x 3.0mm Analytical Column PEEK	
HP-0530-W085	HAISIL 300 C8 5µm, 5cm x 3.0mm Analytical Column PEEK	
HP-0530-W185	HAISIL 300 C18 5µm, 5cm x 3.0mm Analytical Column PEEK	
HP-0546-W045	HAISIL 300 C4 5µm, 5cm x 4.6mm Analytical Column PEEK	
HP-0546-W085	HAISIL 300 C8 5µm, 5cm x 4.6mm Analytical Column PEEK	
HP-0546-W185	HAISIL 300 C18 5µm, 5cm x 4.6mm Analytical Column PEEK	
HP-0721-W045	HAISIL 300 C4 5µm, 7cm x 2.1mm Analytical Column PEEK	
HP-0721-W085	HAISIL 300 C8 5µm, 7cm x 2.1mm Analytical Column PEEK	
HP-0721-W185	HAISIL 300 C18 5µm, 7cm x 2.1mm Analytical Column PEEK	
HP-0730-W045	HAISIL 300 C4 5µm, 7cm x 3.0mm Analytical Column PEEK	
HP-0730-W085	HAISIL 300 C8 5µm, 7cm x 3.0mm Analytical Column PEEK	
HP-0730-W185	HAISIL 300 C18 5µm, 7cm x 3.0mm Analytical Column PEEK	
HP-0746-W045	HAISIL 300 C4 5µm, 7cm x 4.6mm Analytical Column PEEK	
HP-0746-W085	HAISIL 300 C8 5µm, 7cm x 4.6mm Analytical Column PEEK	
HP-0746-W185	HAISIL 300 C18 5µm, 7cm x 4.6mm Analytical Column PEEK	
HP-1021-W045	HAISIL 300 C4 5µm, 10cm x 2.1mm Analytical Column PEEK	
HP-1021-W085	HAISIL 300 C8 5µm, 10cm x 2.1mm Analytical Column PEEK	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
HP-1021-W185	HAISIL 300 C18 5µm, 10cm x 2.1mm Analytical Column PEEK	
HP-1030-W045	HAISIL 300 C4 5µm, 10cm x 3.0mm Analytical Column PEEK	
HP-1030-W085	HAISIL 300 C8 5µm, 10cm x 3.0mm Analytical Column PEEK	
HP-1030-W185	HAISIL 300 C18 5µm, 10cm x 3.0mm Analytical Column PEEK	
HP-1046-W045	HAISIL 300 C4 5µm, 10cm x 4.6mm Analytical Column PEEK	
HP-1046-W085	HAISIL 300 C8 5µm, 10cm x 4.6mm Analytical Column PEEK	
HP-1046-W185	HAISIL 300 C18 5µm, 10cm x 4.6mm Analytical Column PEEK	
HP-1530-W045	HAISIL 300 C4 5µm, 15cm x 3.0mm Analytical Column PEEK	
HP-1530-W085	HAISIL 300 C8 5µm, 15cm x 3.0mm Analytical Column PEEK	
HP-1530-W185	HAISIL 300 C18 5µm, 15cm x 3.0mm Analytical Column PEEK	
HP-1546-W045	HAISIL 300 C4 5µm, 15cm x 4.6mm Analytical Column PEEK	
HP-1546-W085	HAISIL 300 C8 5µm, 15cm x 4.6mm Analytical Column PEEK	
HP-1546-W185	HAISIL 300 C18 5µm, 15cm x 4.6mm Analytical Column PEEK	
HS-0301-M085	HAISIL 100 C8 5µm, 3cm x 1.0mm Microbore Column s.s.	
HS-0301-M183	HAISIL 100 C18 3µm, 3cm x 1.0mm Microbore Column s.s.	
HS-0301-M185	HAISIL 100 C18 5µm, 3cm x 1.0mm Microbore Column s.s.	
HS-0301-SIL5	HAISIL 100 Silica 5µm, 3cm x 1.0mm Microbore Column s.s.	
HS-0301-W045	HAISIL 300 C4 5µm, 3cm x 1.0mm Microbore Column s.s.	
HS-0301-W085	HAISIL 300 C8 5µm, 3cm x 1.0mm Microbore Column s.s.	
HS-0301-W185	HAISIL 300 C18 5µm, 3cm x 1.0mm Microbore Column s.s.	
HS-0321-M085	HAISIL 100 C8 5µm, 3cm x 2.1mm Analytical Column s.s.	
HS-0321-M183	HAISIL 100 C18 3µm, 3cm x 2.1mm Analytical Column s.s.	
HS-0321-M185	HAISIL 100 C18 5µm, 3cm x 2.1mm Analytical Column s.s.	
HS-0321-SIL5	HAISIL 100 Silica 5µm, 3cm x 2.1mm Analytical Column	
HS-0321-W045	HAISIL 300 C4 5µm, 3cm x 2.1mm Analytical Column s.s.	
HS-0321-W085	HAISIL 300 C8 5µm, 3cm x 2.1mm Analytical Column s.s.	
HS-0321-W185	HAISIL 300 C18 5µm, 3cm x 2.1mm Analytical Column s.s.	
HS-0330-M085	HAISIL 100 C8 5µm, 3cm x 3.0mm Analytical Column s.s.	
HS-0330-M183	HAISIL 100 C18 3µm, 3cm x 3.0mm Analytical Column s.s.	
HS-0330-M185	HAISIL 100 C18 5µm, 3cm x 3.0mm Analytical Column s.s.	
HS-0330-SIL5	HAISIL 100 Silica 5µm, 3cm x 3.0mm Analytical Column s.s.	
HS-0330-W045	HAISIL 300 C4 5µm, 3cm x 3.0mm Analytical Column s.s.	
HS-0330-W085	HAISIL 300 C8 5µm, 3cm x 3.0mm Analytical Column s.s.	
HS-0330-W185	HAISIL 300 C18 5µm, 3cm x 3.0mm Analytical Column s.s.	
HS-0346-M085	HAISIL 100 C8 5µm, 3cm x 4.6mm Analytical Column s.s.	
HS-0346-M183	HAISIL 100 C18 3µm, 3cm x 4.6mm Analytical Column s.s.	
HS-0346-M185	HAISIL 100 C18 5µm, 3cm x 4.6mm Analytical Column s.s.	
HS-0346-SIL5	HAISIL 100 Silica 5µm, 3cm x 4.6mm Analytical Column s.s.	
HS-0346-W045	HAISIL 300 C4 5µm, 3cm x 4.6mm Analytical Column s.s.	
HS-0346-W085	HAISIL 300 C8 5µm, 3cm x 4.6mm Analytical Column s.s.	
HS-0346-W185	HAISIL 300 C18 5µm, 3cm x 4.6mm Analytical Column s.s.	
HS-03M3-M085	HAISIL 100 C8 5µm, 3cm x 0.3mm Microbore Column s.s.	
HS-03M3-M183	HAISIL 100 C18 3µm, 3cm x 0.3mm Microbore Column	
HS-03M3-M185	HAISIL 100 C18 5µm, 3cm x 0.3mm Microbore Column s.s.	
HS-03M3-SIL5	HAISIL 100 Silica 5µm, 3cm x 0.3mm Microbore Column	
HS-03M3-W045	HAISIL 300 C4 5µm, 3cm x 0.3mm Microbore Column s.s.	
HS-03M3-W085	HAISIL 300 C8 5µm, 3cm x 0.3mm Microbore Column s.s.	
HS-03M3-W185	HAISIL 300 C18 5µm, 3cm x 0.3mm Microbore Column s.s.	
HS-03M5-M085	HAISIL 100 C8 5µm, 3cm x 0.5mm Microbore Column s.s.	
HS-03M5-M183	HAISIL 100 C18 3µm, 3cm x 0.5mm Microbore Column s.s.	
HS-03M5-M185	HAISIL 100 C18 5µm, 3cm x 0.5mm Microbore Column s.s.	
HS-03M5-SIL5	HAISIL 100 Silica 5µm, 3cm x 0.5mm Microbore Column s.s.	
HS-03M5-W045	HAISIL 300 C4 5µm, 3cm x 0.5mm Microbore Column s.s.	
HS-03M5-W085	HAISIL 300 C8 5µm, 3cm x 0.5mm Microbore Column s.s.	
HS-03M5-W185	HAISIL 300 C18 5µm, 3cm x 0.5mm Microbore Column s.s.	
HS-0501-M085	HAISIL 100 C8 5µm, 5cm x 1.0mm Microbore Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
HS-0501-M183	HAISIL 100 C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
HS-0501-M185	HAISIL 100 C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
HS-0501-SIL5	HAISIL 100 Silica 5µm, 5cm x 1.0mm Microbore Column s.s.	
HS-0501-W045	HAISIL 300 C4 5µm, 5cm x 1.0mm Microbore Column s.s.	
HS-0501-W085	HAISIL 300 C8 5µm, 5cm x 1.0mm Microbore Column s.s.	
HS-0501-W185	HAISIL 300 C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
HS-0510-M085	HAISIL 100 C8 5µm, 5cm x 10mm Semi/Prep Column s.s.	
HS-0510-M185	HAISIL 100 C18 5µm, 5cm x 10mm Semi/Prep Column s.s.	
HS-0510-SIL5	HAISIL 100 Silica 5µm, 5cm x 10mm Semi/Prep Column s.s.	
HS-0510-W045	HAISIL 300 C4 5µm, 5cm x 10mm Semi/Prep Column	
HS-0510-W085	HAISIL 300 C8 5µm, 5cm x 10mm Semi/Prep Column	
HS-0510-W185	HAISIL 300 C18 5µm, 5cm x 10mm Semi/Prep Column	
HS-0520-M085	HAISIL 100 C8 5µm, 5cm x 20mm Semi/Prep Column s.s.	
HS-0520-M185	HAISIL 100 C18 5µm, 5cm x 20mm Semi/Prep Column s.s.	
HS-0520-SIL5	HAISIL 100 Silica 5µm, 5cm x 20mm Semi/Prep Column s.s.	
HS-0520-W045	HAISIL 300 C4 5µm, 5cm x 20mm Semi/Prep Column s.s.	
HS-0520-W085	HAISIL 300 C8 5µm, 5cm x 20mm Semi/Prep Column s.s.	
HS-0520-W185	HAISIL 300 C18 5µm, 5cm x 20mm Semi/Prep Column s.s.	
HS-0521-M085	HAISIL 100 C8 5µm, 5cm x 2.1mm Analytical Column s.s.	
HS-0521-M183	HAISIL 100 C18 3µm, 5cm x 2.1mm Analytical Column s.s.	
HS-0521-M185	HAISIL 100 C18 5µm, 5cm x 2.1mm Analytical Column s.s.	
HS-0521-SIL5	HAISIL 100 Silica 5µm, 5cm x 2.1mm Analytical Column s.s.	
HS-0521-W045	HAISIL 300 C4 5µm, 5cm x 2.1mm Analytical Column s.s.	
HS-0521-W085	HAISIL 300 C8 5µm, 5cm x 2.1mm Analytical Column s.s.	
HS-0521-W185	HAISIL 300 C18 5µm, 5cm x 2.1mm Analytical Column s.s.	
HS-0530-M085	HAISIL 100 C8 5µm, 5cm x 3.0mm Analytical Column s.s.	
HS-0530-M183	HAISIL 100 C18 3µm, 5cm x 3.0mm Analytical Column s.s.	
HS-0530-M185	HAISIL 100 C18 5µm, 5cm x 3.0mm Analytical Column s.s.	
HS-0530-SIL5	HAISIL 100 Silica 5µm, 5cm x 3.0mm Analytical Column s.s.	
HS-0530-W045	HAISIL 300 C4 5µm, 5cm x 3.0mm Analytical Column s.s.	
HS-0530-W085	HAISIL 300 C8 5µm, 5cm x 3.0mm Analytical Column s.s.	
HS-0530-W185	HAISIL 300 C18 5µm, 5cm x 3.0mm Analytical Column s.s.	
HS-0546-M085	HAISIL 100 C8 5µm, 5cm x 4.6mm Analytical Column s.s.	
HS-0546-M183	HAISIL 100 C18 3µm, 5cm x 4.6mm Analytical Column s.s.	
HS-0546-M185	HAISIL 100 C18 5µm, 5cm x 4.6mm Analytical Column s.s.	
HS-0546-SIL5	HAISIL 100 Silica 5µm, 5cm x 4.6mm Analytical Column s.s.	
HS-0546-W045	HAISIL 300 C4 5µm, 5cm x 4.6mm Analytical Column s.s.	
HS-0546-W085	HAISIL 300 C8 5µm, 5cm x 4.6mm Analytical Column s.s.	
HS-0546-W185	HAISIL 300 C18 5µm, 5cm x 4.6mm Analytical Column s.s.	
HS-05M3-M085	HAISIL 100 C8 5µm, 5cm x 0.3mm Microbore Column s.s.	
HS-05M3-M183	HAISIL 100 C18 3µm, 5cm x 0.3mm Microbore Column s.s.	
HS-05M3-M185	HAISIL 100 C18 5µm, 5cm x 0.3mm Microbore Column s.s.	
HS-05M3-SIL5	HAISIL 100 Silica 5µm, 5cm x 0.3mm Microbore Column	
HS-05M3-W045	HAISIL 300 C4 5µm, 5cm x 0.3mm Microbore Column s.s.	
HS-05M3-W085	HAISIL 300 C8 5µm, 5cm x 0.3mm Microbore Column s.s.	
HS-05M3-W185	HAISIL 300 C18 5µm, 5cm x 0.3mm Microbore Column s.s.	
HS-05M5-M085	HAISIL 100 C8 5µm, 5cm x 0.5mm Microbore Column s.s.	
HS-05M5-M183	HAISIL 100 C18 3µm, 5cm x 0.5mm Microbore Column s.s.	
HS-05M5-M185	HAISIL 100 C18 5µm, 5cm x 0.5mm Microbore Column s.s.	
HS-05M5-SIL5	HAISIL 100 Silica 5µm, 5cm x 0.5mm Microbore Column s.s.	
HS-05M5-W045	HAISIL 300 C4 5µm, 5cm x 0.5mm Microbore Column s.s.	
HS-05M5-W085	HAISIL 300 C8 5µm, 5cm x 0.5mm Microbore Column s.s.	
HS-05M5-W185	HAISIL 300 C18 5µm, 5cm x 0.5mm Microbore Column s.s.	
HS-1001-M085	HAISIL 100 C8 5µm, 10cm x 1.0mm Microbore Column s.s.	
HS-1001-M183	HAISIL 100 C18 3µm, 10cm x 1.0mm Microbore Column s.s.	
HS-1001-M185	HAISIL 100 C18 5µm, 10cm x 1.0mm Microbore Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
HS-1001-SIL5	HAISIL 100 Silica 5µm, 10cm x 1.0mm Microbore Column s.s.	
HS-1001-W045	HAISIL 300 C4 5µm, 10cm x 1.0mm Microbore Column s.s.	
HS-1001-W085	HAISIL 300 C8 5µm, 10cm x 1.0mm Microbore Column s.s.	
HS-1001-W185	HAISIL 300 C18 5µm, 10cm x 1.0mm Microbore Column s.s.	
HS-1010-M085	HAISIL 100 C8 5µm, 10cm x 10mm Semi/Prep Column s.s.	
HS-1010-M185	HAISIL 100 C18 5µm, 10cm x 10mm Semi/Prep Column s.s.	
HS-1010-SIL5	HAISIL 100 Silica 5µm, 10cm x 10mm Semi/Prep Column s.s.	
HS-1010-W045	HAISIL 300 C4 5µm, 10cm x 10mm Semi/Prep Column s.s.	
HS-1010-W085	HAISIL 300 C8 5µm, 10cm x 10mm Semi/Prep Column s.s.	
HS-1010-W185	HAISIL 300 C18 5µm, 10cm x 10mm Semi/Prep Column s.s.	
HS-1020-M085	HAISIL 100 C8 5µm, 10cm x 20mm Semi/Prep Column s.s.	
HS-1020-M185	HAISIL 100 C18 5µm, 10cm x 20mm Semi/Prep Column s.s.	
HS-1020-SIL5	HAISIL 100 Silica 5µm, 10cm x 20mm Semi/Prep Column s.s.	
HS-1020-W045	HAISIL 300 C4 5µm, 10cm x 20mm Semi/Prep Column s.s.	
HS-1020-W085	HAISIL 300 C8 5µm, 10cm x 20mm Semi/Prep Column s.s.	
HS-1020-W185	HAISIL 300 C18 5µm, 10cm x 20mm Semi/Prep Column s.s.	
HS-1021-M085	HAISIL 100 C8 5µm, 10cm x 2.1mm Analytical Column s.s.	
HS-1021-M183	HAISIL 100 C18 3µm, 10cm x 2.1mm Analytical Column s.s.	
HS-1021-M185	HAISIL 100 C18 5µm, 10cm x 2.1mm Analytical Column s.s.	
HS-1021-SIL5	HAISIL 100 Silica 5µm, 10cm x 2.1mm Analytical Column s.s.	
HS-1021-W045	HAISIL 300 C4 5µm, 10cm x 2.1mm Analytical Column s.s.	
HS-1021-W085	HAISIL 300 C8 5µm, 10cm x 2.1mm Analytical Column s.s.	
HS-1021-W185	HAISIL 300 C18 5µm, 10cm x 2.1mm Analytical Column s.s.	
HS-1030-M085	HAISIL 100 C8 5µm, 10cm x 3.0mm Analytical Column s.s.	
HS-1030-M183	HAISIL 100 C18 3µm, 10cm x 3.0mm Analytical Column s.s.	
HS-1030-M185	HAISIL 100 C18 5µm, 10cm x 3.0mm Analytical Column s.s.	
HS-1030-SIL5	HAISIL 100 Silica 5µm, 10cm x 3.0mm Analytical Column s.s.	
HS-1030-W045	HAISIL 300 C4 5µm, 10cm x 3.0mm Analytical Column s.s.	
HS-1030-W085	HAISIL 300 C8 5µm, 10cm x 3.0mm Analytical Column s.s.	
HS-1030-W185	HAISIL 300 C18 5µm, 10cm x 3.0mm Analytical Column s.s.	
HS-1046-M085	HAISIL 100 C8 5µm, 10cm x 4.6mm Analytical Column s.s.	
HS-1046-M183	HAISIL 100 C18 3µm, 10cm x 4.6mm Analytical Column s.s.	
HS-1046-M185	HAISIL 100 C18 5µm, 10cm x 4.6mm Analytical Column s.s.	
HS-1046-SIL5	HAISIL 100 Silica 5µm, 10cm x 4.6mm Analytical Column s.s.	
HS-1046-W045	HAISIL 300 C4 5µm, 10cm x 4.6mm Analytical Column s.s.	
HS-1046-W085	HAISIL 300 C8 10cm x 4.6mm Analytical Column s.s.	
HS-1046-W185	HAISIL 300 C18 5µm, 10cm x 4.6mm Analytical Column s.s.	
HS-10M3-M085	HAISIL 100 C8 5µm, 10cm x 0.3mm Microbore Column s.s.	
HS-10M3-M183	HAISIL 100 C18 3µm, 10cm x 0.3mm Microbore Column s.s.	
HS-10M3-M185	HAISIL 100 C18 5µm, 10cm x 0.3mm Microbore Column s.s.	
HS-10M3-SIL5	HAISIL 100 Silica 5µm, 10cm x 0.3mm Microbore Column s.s.	
HS-10M3-W045	HAISIL 300 C4 5µm, 10cm x 0.3mm Microbore Column s.s.	
HS-10M3-W085	HAISIL 300 C8 5µm, 10cm x 0.3mm Microbore Column s.s.	
HS-10M3-W185	HAISIL 300 C18 5µm, 10cm x 0.3mm Microbore Column s.s.	
HS-10M5-M085	HAISIL 100 C8 5µm, 10cm x 0.5mm Microbore Column s.s.	
HS-10M5-M183	HAISIL 100 C18 3µm, 10cm x 0.5mm Microbore Column s.s.	
HS-10M5-M185	HAISIL 100 C18 5µm, 10cm x 0.5mm Microbore Column s.s.	
HS-10M5-SIL5	HAISIL 100 Silica 5µm, 10cm x 0.5mm Microbore Column s.s.	
HS-10M5-W045	HAISIL 300 C4 5µm, 10cm x 0.5mm Microbore Column s.s.	
HS-10M5-W085	HAISIL 300 C8 5µm, 10cm x 0.5mm Microbore Column s.s.	
HS-10M5-W185	HAISIL 300 C18 5µm, 10cm x 0.5mm Microbore Column s.s.	
HS-1501-M085	HAISIL 100 C8 5µm, 15cm x 1.0mm Microbore Column s.s.	
HS-1501-M183	HAISIL 100 C18 3µm, 15cm x 1.0mm Microbore Column s.s.	
HS-1501-M185	HAISIL 100 C18 5µm, 15cm x 1.0mm Microbore Column s.s.	
HS-1501-SIL5	HAISIL 100 Silica 5µm, 15cm x 1.0mm Microbore Column s.s.	
HS-1501-W045	HAISIL 300 C4 5µm, 15cm x 1.0mm Microbore Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
HS-1501-W085	HAISIL 300 C8 5µm, 15cm x 1.0mm Microbore Column s.s.	
HS-1501-W185	HAISIL 300 C18 5µm, 15cm x 1.0mm Microbore Column s.s.	
HS-1510-M085	HAISIL 100 C8 5µm, 15cm x 10mm Semi/Prep Column s.s.	
HS-1510-M185	HAISIL 100 C18 5µm, 15cm x 10mm Semi/Prep Column s.s.	
HS-1510-SIL5	HAISIL 100 Silica 5µm, 15cm x 10mm Semi/Prep Column s.s.	
HS-1510-W045	HAISIL 300 C4 5µm, 15cm x 10mm Semi/Prep Column s.s.	
HS-1510-W085	HAISIL 300 C8 5µm, 15cm x 10mm Semi/Prep Column s.s.	
HS-1510-W185	HAISIL 300 C18 5µm, 15cm x 10mm Semi/Prep Column s.s.	
HS-1520-M085	HAISIL 100 C8 5µm, 15cm x 20mm Semi/Prep Column s.s.	
HS-1520-M185	HAISIL 100 C18 5µm, 15cm x 20mm Semi/Prep Column s.s.	
HS-1520-SIL5	HAISIL 100 Silica 5µm, 15cm x 20mm Semi/Prep Column s.s.	
HS-1520-W045	HAISIL 300 C4 5µm, 15cm x 20mm Semi/Prep Column s.s.	
HS-1520-W085	HAISIL 300 C8 5µm, 15cm x 20mm Semi/Prep Column s.s.	
HS-1520-W185	HAISIL 300 C18 5µm, 15cm x 20mm Semi/Prep Column s.s.	
HS-1521-M085	HAISIL 100 C8 5µm, 15cm x 2.1mm Analytical Column s.s.	
HS-1521-M183	HAISIL 100 C18 3µm, 15cm x 2.1mm Analytical Column s.s.	
HS-1521-M185	HAISIL 100 C18 5µm, 15cm x 2.1mm Analytical Column s.s.	
HS-1521-SIL5	HAISIL 100 Silica 5µm, 15cm x 2.1mm Analytical Column s.s.	
HS-1521-W045	HAISIL 300 C4 5µm, 15cm x 2.1mm Analytical Column s.s.	
HS-1521-W085	HAISIL 300 C8 5µm, 15cm x 2.1mm Analytical Column s.s.	
HS-1521-W185	HAISIL 300 C18 5µm, 15cm x 2.1mm Analytical Column s.s.	
HS-1530-M085	HAISIL 100 C8 5µm, 15cm x 3.0mm Analytical Column s.s.	
HS-1530-M183	HAISIL 100 C18 3µm, 15cm x 3.0mm Analytical Column s.s.	
HS-1530-M185	HAISIL 100 C18 5µm, 15cm x 3.0mm Analytical Column s.s.	
HS-1530-SIL5	HAISIL 100 Silica 5µm, 15cm x 3.0mm Analytical Column s.s.	
HS-1530-W045	HAISIL 300 C4 5µm, 15cm x 3.0mm Analytical Column s.s.	
HS-1530-W085	HAISIL 300 C8 5µm, 15cm x 3.0mm Analytical Column s.s.	
HS-1530-W185	HAISIL 300 C18 5µm, 15cm x 3.0mm Analytical Column s.s.	
HS-1546-M085	HAISIL 100 C8 5µm, 15cm x 4.6mm Analytical Column s.s.	
HS-1546-M183	HAISIL 100 C18 3µm, 15cm x 4.6mm Analytical Column s.s.	
HS-1546-M185	HAISIL 100 C18 5µm, 15cm x 4.6mm Analytical Column s.s.	
HS-1546-SIL5	HAISIL 100 Silica 5µm, 15cm x 4.6mm Analytical Column s.s.	
HS-1546-W045	HAISIL 300 C4 5µm, 15cm x 4.6mm Analytical Column s.s.	
HS-1546-W085	HAISIL 300 C8 5µm, 15cm x 4.6mm Analytical Column s.s.	
HS-1546-W185	HAISIL 300 C18 5µm, 15cm x 4.6mm Analytical Column s.s.	
HS-15M5-M085	HAISIL 100 C8 5µm, 15cm x 0.5mm Microbore Column s.s.	
HS-15M5-M183	HAISIL 100 C18 3µm, 15cm x 0.5mm Microbore Column s.s.	
HS-15M5-M185	HAISIL 100 C18 5µm, 15cm x 0.5mm Microbore Column s.s.	
HS-15M5-SIL5	HAISIL 100 Silica 5µm, 15cm x 0.5mm Microbore Column s.s.	
HS-15M5-W045	HAISIL 300 C4 5µm, 15cm x 0.5mm Microbore Column s.s.	
HS-15M5-W085	HAISIL 300 C8 5µm, 15cm x 0.5mm Microbore Column s.s.	
HS-15M5-W185	HAISIL 300 C18 5µm, 15cm x 0.5mm Microbore Column s.s.	
HS-2501-M085	HAISIL 100 C8 5µm, 25cm x 1.0mm Microbore Column s.s.	
HS-2501-M183	HAISIL 100 C18 3µm, 25cm x 1.0mm Microbore Column s.s.	
HS-2501-M185	HAISIL 100 C18 5µm, 25cm x 1.0mm Microbore Column s.s.	
HS-2501-SIL5	HAISIL 100 Silica 5µm, 25cm x 1.0mm Microbore Column s.s.	
HS-2501-W045	HAISIL 300 C4 5µm, 25cm x 1.0mm Microbore Column s.s.	
HS-2501-W085	HAISIL 300 C8 5µm, 25cm x 1.0mm Microbore Column s.s.	
HS-2501-W185	HAISIL 300 C18 5µm, 25cm x 1.0mm Microbore Column s.s.	
HS-2510-M085	HAISIL 100 C8 5µm, 25cm x 10mm Semi/Prep Column s.s.	
HS-2510-M185	HAISIL 100 C18 5µm, 25cm x 10mm Semi/Prep Column s.s.	
HS-2510-SIL5	HAISIL 100 Silica 5µm, 25cm x 10mm Semi/Prep Column s.s.	
HS-2510-W045	HAISIL 300 C4 5µm, 25cm x 10mm Semi/Prep Column s.s.	
HS-2510-W085	HAISIL 300 C8 5µm, 25cm x 10mm Semi/Prep Column s.s.	
HS-2510-W185	HAISIL 300 C18 5µm, 25cm x 10mm Semi/Prep Column s.s.	
HS-2520-M085	HAISIL 100 C8 5µm, 25cm x 20mm Semi/Prep Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
HS-2520-M185	HAISIL 100 C18 5µm, 25cm x 20mm Semi/Prep Column s.s.	
HS-2520-SIL5	HAISIL 100 Silica 5µm, 25cm x 20mm Semi/Prep Column s.s.	
HS-2520-W045	HAISIL 300 C4 5µm, 25cm x 20mm Semi/Prep Column s.s.	
HS-2520-W085	HAISIL 300 C8 5µm, 25cm x 20mm Semi/Prep Column s.s.	
HS-2520-W185	HAISIL 300 C18 5µm, 25cm x 20mm Semi/Prep Column s.s.	
HS-2521-M085	HAISIL 100 C8 5µm, 25cm x 2.1mm Analytical Column s.s.	
HS-2521-M183	HAISIL 100 C18 3µm, 25cm x 2.1mm Analytical Column s.s.	
HS-2521-M185	HAISIL 100 C18 5µm, 25cm x 2.1mm Analytical Column s.s.	
HS-2521-SIL5	HAISIL 100 Silica 5µm, 25cm x 2.1mm Analytical Column s.s.	
HS-2521-W045	HAISIL 300 C4 5µm, 25cm x 2.1mm Analytical Column s.s.	
HS-2521-W085	HAISIL 300 C8 5µm, 25cm x 2.1mm Analytical Column s.s.	
HS-2521-W185	HAISIL 300 C18 5µm, 25cm x 2.1mm Analytical Column s.s.	
HS-2530-M085	HAISIL 100 C8 5µm, 25cm x 3.0mm Analytical Column s.s.	
HS-2530-M183	HAISIL 100 C18 3µm, 25cm x 3.0mm Analytical Column s.s.	
HS-2530-M185	HAISIL 100 C18 5µm, 25cm x 3.0mm Analytical Column s.s.	
HS-2530-SIL5	HAISIL 100 Silica 5µm, 25cm x 3.0mm Analytical Column s.s.	
HS-2530-W045	HAISIL 300 C4 5µm, 25cm x 3.0mm Analytical Column s.s.	
HS-2530-W085	HAISIL 300 C8 5µm, 25cm x 3.0mm Analytical Column s.s.	
HS-2530-W185	HAISIL 300 C18 5µm, 25cm x 3.0mm Analytical Column s.s.	
HS-2546-M085	HAISIL 100 C8 5µm, 25cm x 4.6mm Analytical Column s.s.	
HS-2546-M183	HAISIL 100 C18 3µm, 25cm x 4.6mm Analytical Column s.s.	
HS-2546-M185	HAISIL 100 C18 5µm, 25cm x 4.6mm Analytical Column s.s.	
HS-2546-SIL5	HAISIL 100 Silica 5µm, 25cm x 4.6mm Analytical Column s.s.	
HS-2546-W045	HAISIL 300 C4 5µm, 25cm x 4.6mm Analytical Column s.s.	
HS-2546-W085	HAISIL 300 C8 5µm, 25cm x 4.6mm Analytical Column s.s.	
HS-2546-W185	HAISIL 300 C18 5µm, 25cm x 4.6mm Analytical Column s.s.	
HS-25M5-M085	HAISIL 100 C8 5µm, 25cm x 0.5mm Microbore Column s.s.	
HS-25M5-M183	HAISIL 100 C18 3µm, 25cm x 0.5mm Microbore Column s.s.	
HS-25M5-M185	HAISIL 100 C18 5µm, 25cm x 0.5mm Microbore Column s.s.	
HS-25M5-SIL5	HAISIL 100 Silica 5µm, 25cm x 0.5mm Microbore Column s.s.	
K108	Kromasil C8 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
K118	Kromasil C18 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
KF-0310-C085	Kromasil C8 5µm, 3cm x 10mm HAlspeed Cartridge s.s. 1/pkg	
KF-0310-C185	Kromasil C18 5µm, 3cm x 10mm HAlspeed Cartridge s.s. 1/pkg	
KF-0432-C085	Kromasil C8 5µm, 4cm x 3.2mm Analytical Cartridge s.s. 2/pkg	
KF-0432-C183	Kromasil C18 3µm, 4cm x 3.2mm Analytical Cartridge 2/pkg	
KF-0432-C185	Kromasil C18 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
KF-0446-C085	Kromasil C8 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
KF-0446-C183	Kromasil C18 3µm, 4cm x 4.6mm Analytical Cartridge 2/pkg	
KF-0446-C185	Kromasil C18 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
KF-0510-C085	Kromasil C8 5µm, 5cm x 10mm HAlspeed Cartridge 1/pkg	
KF-0510-C185	Kromasil C18 5µm, 5cm x 10mm HAlspeed Cartridge 1/pkg	
KF-0532-C085	Kromasil C8 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
KF-0532-C183	Kromasil C18 3µm, 5cm x 3.2mm Analytical Cartridge 2/pkg	
KF-0532-C185	Kromasil C18 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
KF-0546-C085	Kromasil C8 5µm, 5cm x 4.6mm Analytical Cartridge s.s. 2/pkg	
KF-0546-C183	Kromasil C18 3µm, 5cm 4.6mm Analytical Cartridge 2/pkg	
KF-0546-C185	Kromasil C18 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
KF-1010-C085	Kromasil C8 5µm, 10cm x 10mm HAlspeed Cartridge 1/pkg	
KF-1010-C185	Kromasil C18 5µm, 10cm x 10mm HAlSpeed Cartridge 1/pkg	
KF-1032-C085	Kromasil C8 5µm, 10cm x 3.2mm Cartridge s.s.	
KF-1032-C183	Kromasil C18 3µm, 10cm x3.2mm Analytical Cartridge 1/pkg	
KF-1032-C185	Kromasil C18 5µm, 10cm x 3.2mm Cartridge s.s.	
KF-1046-C085	Kromasil C8 5µm, 10cm x 4.6mm Cartridge s.s.	
KF-1046-C183	Kromasil C18 3µm, 10cm x 4.6mm Analytical Cartridge 1/pkg	
KF-1046-C185	Kromasil C18 5µm, 10cm x 4.6mm Cartridge s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
KK-0221-C085	Kromasil C8 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
KK-0221-C185	Kromasil C18 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
KK-0421-C085	Kromasil C8 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
KK-0421-C185	Kromasil C18 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
KK-0521-C085	Kromasil C8 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
KK-0521-C185	Kromasil C18 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
KS-0301-C085	Kromasil C8 5µm, 3cm x 1.0mm Microbore Column s.s.	
KS-0301-C183	Kromasil C18 3µm, 3cm x 1.0mm Microbore Column s.s.	
KS-0301-C185	Kromasil C18 5µm, 3cm x 1.0mm Microbore Column s.s.	
KS-0321-C085	Kromasil C8 5µm, 3cm x 2.1mm Analytical Column s.s.	
KS-0321-C183	Kromasil C18 3µm, 3cm x 2.1mm Analytical Column s.s.	
KS-0321-C185	Kromasil C18 5µm, 3cm x 2.1mm Analytical Column s.s.	
KS-0330-C085	Kromasil C8 5µm, 3cm x 3.0mm Analytical Column s.s.	
KS-0330-C183	Kromasil C18 3µm, 3cm x 3.0mm Analytical Column s.s.	
KS-0330-C185	Kromasil C18 5µm, 3cm x 3.0mm Analytical Column s.s.	
KS-0346-C085	Kromasil C8 5µm, 3cm x 4.6mm Analytical Column s.s.	
KS-0346-C183	Kromasil C18 3µm, 3cm x 4.6mm Analytical Column s.s.	
KS-0346-C185	Kromasil C18 5µm, 3cm x 4.6mm Analytical Column s.s.	
KS-03M3-C085	Kromasil C8 5µm, 3cm x 0.3mm Microbore Column s.s.	
KS-03M3-C183	Kromasil C18 3µm, 3cm x 0.3mm Microbore Column s.s.	
KS-03M3-C185	Kromasil C18 5µm, 3cm x 0.3mm Microbore Column s.s.	
KS-03M5-C085	Kromasil C8 5µm, 3cm x 0.5mm Microbore Column s.s.	
KS-03M5-C183	Kromasil C18 3µm, 3cm x 0.5mm Microbore Column s.s.	
KS-03M5-C185	Kromasil C18 5µm, 3cm x 0.5mm Microbore Column s.s.	
KS-0501-C085	Kromasil C8 5µm, 5cm x 1.0mm Microbore Column s.s.	
KS-0501-C183	Kromasil C18 3µm, 5cm x 1.0mm Microbore Column s.s.	
KS-0501-C185	Kromasil C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
KS-0510-C085	Kromasil C8 5µm, 5cm x 10mm Semi/Prep Column s.s.	
KS-0510-C185	Kromasil C18 5µm, 5cm x 10mm Semi/Prep Column s.s.	
KS-0520-C085	Kromasil C8 5µm, 5cm x 20mm Semi/Prep Column s.s.	
KS-0520-C185	Kromasil C18 5µm, 5cm x 20mm Semi/Prep Column s.s.	
KS-0521-C085	Kromasil C8 5µm, 5cm x 2.1mm Analytical Column s.s.	
KS-0521-C183	Kromasil C18 3µm, 5cm x 2.1mm Analytical Column s.s.	
KS-0521-C185	Kromasil C18 5µm, 5cm x 2.1mm Analytical Column s.s.	
KS-0530-C085	Kromasil C8 5µm, 5cm x 3.0mm Analytical Column s.s.	
KS-0530-C183	Kromasil C18 3µm, 5cm x 3.0mm Analytical Column s.s.	
KS-0530-C185	Kromasil C18 5µm, 5cm x 3.0mm Analytical Column s.s.	
KS-0546-C085	Kromasil C8 5µm, 5cm x 4.6mm Analytical Column s.s.	
KS-0546-C183	Kromasil C18 3µm, 5cm x 4.6mm Analytical Column s.s.	
KS-0546-C185	Kromasil C18 5µm, 5cm x 4.6mm Analytical Column s.s.	
KS-05M3-C085	Kromasil C8 5µm, 5cm x 0.3mm Microbore Column s.s.	
KS-05M3-C183	Kromasil C18 3µm, 5cm x 0.3mm Microbore Column s.s.	
KS-05M3-C185	Kromasil C18 5µm, 5cm x 0.3mm Microbore Column s.s.	
KS-05M5-C085	Kromasil C8 5µm, 5cm x 0.5mm Microbore Column s.s.	
KS-05M5-C183	Kromasil C18 3µm, 5cm x 0.5mm Microbore Column s.s.	
KS-05M5-C185	Kromasil C18 5µm, 5cm x 0.5mm Microbore Column s.s.	
KS-1001-C085	Kromasil C8 5µm, 10cm x 1.0mm Microbore Column s.s.	
KS-1001-C183	Kromasil C18 3µm, 10cm x 1.0mm Microbore Column s.s.	
KS-1001-C185	Kromasil C18 5µm, 10cm x 1.0mm Microbore Column s.s.	
KS-1010-C085	Kromasil C8 5µm, 10cm x 10mm Semi/Prep Column s.s.	
KS-1010-C185	Kromasil C18 5µm, 10cm x 10mm Semi/Prep Column s.s.	
KS-1020-C085	Kromasil C8 5µm, 10cm x 20mm Semi/Prep Column s.s.	
KS-1020-C185	Kromasil C18 5µm, 10cm x 20mm Semi/Prep Column s.s.	
KS-1021-C085	Kromasil C8 5µm, 10cm x 2.1mm Analytical Column s.s.	
KS-1021-C183	Kromasil C18 3µm, 10cm x 2.1mm Analytical Column s.s.	
KS-1021-C185	Kromasil C18 5µm, 10cm x 2.1mm Analytical Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
KS-1030-C085	Kromasil C8 5µm, 10cm x 3.0mm Analytical Column s.s.	
KS-1030-C183	Kromasil C18 3µm, 10cm x 3.0mm Analytical Column s.s.	
KS-1030-C185	Kromasil C18 5µm, 10cm x 3.0mm Analytical Column s.s.	
KS-1046-C085	Kromasil C8 5µm, 10cm x 4.6mm Analytical Column s.s.	
KS-1046-C183	Kromasil C18 3µm, 10cm x 4.6mm Analytical Column s.s.	
KS-1046-C185	Kromasil C18 5µm, 10cm x 4.6mm Analytical Column s.s.	
KS-10M3-C085	Kromasil C8 5µm, 10cm x 0.3mm Microbore Column s.s.	
KS-10M3-C183	Kromasil C18 3µm, 10cm x 0.3mm Microbore Column s.s.	
KS-10M3-C185	Kromasil C18 5µm, 10cm x 0.3mm Microbore Column s.s.	
KS-10M5-C085	Kromasil C8 5µm, 10cm x 0.5mm Microbore Column s.s.	
KS-10M5-C183	Kromasil C18 3µm, 10cm x 0.5mm Microbore Column s.s.	
KS-10M5-C185	Kromasil C18 5µm, 10cm x 0.5mm Microbore Column s.s.	
KS-1501-C085	Kromasil C8 5µm, 15cm x 1.0mm Microbore Column s.s.	
KS-1501-C183	Kromasil C18 3µm, 15cm x 1.0mm Microbore Column s.s.	
KS-1501-C185	Kromasil C18 5µm, 15cm x 1.0mm Microbore Column s.s.	
KS-1510-C085	Kromasil C8 5µm, 15cm x 10mm Semi/Prep Column s.s.	
KS-1510-C185	Kromasil C18 5µm, 15cm x 10mm Semi/Prep Column s.s.	
KS-1520-C085	Kromasil C8 5µm, 15cm x 20mm Semi/Prep Column s.s.	
KS-1520-C185	Kromasil C18 5µm, 15cm x 20mm Semi/Prep Column s.s.	
KS-1521-C085	Kromasil C8 5µm, 15cm x 2.1mm Analytical Column s.s.	
KS-1521-C183	Kromasil C18 3µm, 15cm x 2.1mm Analytical Column s.s.	
KS-1521-C185	Kromasil C18 5µm, 15cm x 2.1mm Analytical Column s.s.	
KS-1530-C085	Kromasil C8 5µm, 15cm x 3.0mm Analytical Column s.s.	
KS-1530-C183	Kromasil C18 3µm, 15cm x 3.0mm Analytical Column s.s.	
KS-1530-C185	Kromasil C18 5µm, 15cm x 3.0mm Analytical Column s.s.	
KS-1546-C085	Kromasil C8 5µm, 15cm x 4.6mm Analytical Column s.s.	
KS-1546-C183	Kromasil C18 3µm, 15cm x 4.6mm Analytical Column s.s.	
KS-1546-C185	Kromasil C18 5µm, 15cm x 4.6mm Analytical Column s.s.	
KS-15M5-C085	Kromasil C8 5µm, 15cm x 0.5mm Microbore Column s.s.	
KS-15M5-C183	Kromasil C18 3µm, 15cm x 0.5mm Microbore Column s.s.	
KS-15M5-C185	Kromasil C18 5µm, 15cm x 0.5mm Microbore Column s.s.	
KS-2501-C085	Kromasil C8 5µm, 25cm x 1.0mm Microbore Column s.s.	
KS-2501-C183	Kromasil C18 3µm, 25cm x 1.0mm Microbore Column s.s.	
KS-2501-C185	Kromasil C18 5µm, 25cm x 1.0mm Microbore Column s.s.	
KS-2510-C085	Kromasil C8 5µm, 25cm x 10mm Semi/Prep Column s.s.	
KS-2510-C185	Kromasil C18 5µm, 25cm x 10mm Semi/Prep Column s.s.	
KS-2520-C085	Kromasil C8 5µm, 25cm x 20mm Semi/Prep Column s.s.	
KS-2520-C185	Kromasil C18 5µm, 25cm x 20mm Semi/Prep Column s.s.	
KS-2521-C085	Kromasil C8 5µm, 25cm x 2.1mm Analytical Column s.s.	
KS-2521-C183	Kromasil C18 3µm, 25cm x 2.1mm Analytical Column s.s.	
KS-2521-C185	Kromasil C18 5µm, 25cm x 2.1mm Analytical Column s.s.	
KS-2530-C085	Kromasil C8 5µm, 25cm x 3.0mm Analytical Column s.s.	
KS-2530-C183	Kromasil C18 3µm, 25cm x 3.0mm Analytical Column s.s.	
KS-2530-C185	Kromasil C18 5µm, 25cm x 3.0mm Analytical Column s.s.	
KS-2546-C085	Kromasil C8 5µm, 25cm x 4.6mm Analytical Column s.s.	
KS-2546-C183	Kromasil C18 3µm, 25cm x 4.6mm Analytical Column s.s.	
KS-2546-C185	Kromasil C18 5µm, 25cm x 4.6mm Analytical Column s.s.	
LI18	HAISIL HL C18 5µm, 2cm x 3.2mm Cartridge s.s. 5/pkg	
LC-0507-C185	HAISIL HL Capellini C18 5µm, 50 x 0.75mm column	
LC-0515-C185	HAISIL HL Capellini C18 5µm, 50 x 0.150mm column	
LC-1007-C185	HAISIL HL Capellini C18 5µm, 100 x 0.075mm column	
LC-1015-C185	HAISIL HL Capellini C18 5µm, 100 x 0.150mm column	
LC-1507-C185	HAISIL HL Capellini C18 5µm, 150 x 0.075mm column	
LC-1515-C185	HAISIL HL Capellini C18 5µm, 150 x 0.150mm column	
LC-2507-C185	HAISIL HL Capellini C18 5µm, 250 x 0.075mm column	
LC-2515-C185	HAISIL HL Capellini C18 5µm, 250 x 0.150mm column	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
LF-0310-C185	HAISIL HL C18 5µm, 3cm x 10mm Cartridge s.s. 1/pkg	
LF-0432-C183	HAISIL HL C18 3µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
LF-0432-C185	HAISIL HL C18 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
LF-0446-C183	HAISIL HL C18 3µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
LF-0446-C185	HAISIL HL C18 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
LF-0510-C185	HAISIL HL C18 5µm, 5cm x 10mm Cartridge s.s. 1/pkg	
LF-0532-C183	HAISIL HL C18 3µm, 5cm 3.2mm Cartridge s.s. 2/pkg	
LF-0532-C185	HAISIL HL C18 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
LF-0546-C183	HAISIL HL C18 3µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
LF-0546-C185	HAISIL HL C18 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
LF-1010-C185	HAISIL HL C18 5µm, 10cm x 10mm Cartridge s.s.	
LF-1032-C183	HAISIL HL C18 3µm, 10cm x 3.2mm Cartridge s.s. 1/pkg	
LF-1032-C185	HAISIL HL C18 5µm, 10cm x 3.2mm Cartridge s.s.	
LF-1046-C183	HAISIL HL C18 3µm, 10cm x 4.6mm Cartridge s.s. 1/pkg	
LF-1046-C185	HAISIL HL C18 5µm, 10cm x 4.6mm Cartridge s.s.	
LK-0221-C185	HAISIL HL C18 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
LK-0421-C185	HAISIL HL C18 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
LK-0521-C185	HAISIL HL C18 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
LP-0101-C185	HAISIL HL Piccolo C18 5µm, 10 x 1.0mm column	
LP-M2M5-C185	HAISIL HL Piccolo C18 5µm, 2.5 x 0.5mm column	
LP-M501-C185	HAISIL HL Piccolo C18 5µm, 5 x 1.0mm column	
LR-0121-C185	HAISIL HL Sprite C18 5µm, 10 x 2.1mm column	
LR-0221-C185	HAISIL HL Sprite C18 5µm, 20 x 2.1mm column	
LR-0421-C185	HAISIL HL Sprite C18 5µm, 40 x 2.1mm column	
LS-0301-C183	HAISIL HL C18 3µm, 3cm x 1.0mm Microbore Column s.s.	
LS-0301-C185	HAISIL HL C18 5µm, 3cm x 1.0mm Microbore Column s.s.	
LS-0321-C183	HAISIL HL C18 3µm, 3cm x 2.1mm Analytical Column s.s.	
LS-0321-C185	HAISIL HL C18 5µm, 3cm x 2.1mm Analytical Column s.s.	
LS-0330-C183	HAISIL HL C18 3µm, 3cm x 3.0mm Analytical Column s.s.	
LS-0330-C185	HAISIL HL C18 5µm, 3cm x 3.0mm Analytical Column s.s.	
LS-0346-C183	HAISIL HL C18 3µm, 3cm x 4.6mm Analytical Column s.s.	
LS-0346-C185	HAISIL HL C18 5µm, 3cm x 4.6mm Analytical Column s.s.	
LS-03M3-C183	HAISIL HL C18 3µm, 3cm x 0.3mm Microbore Column s.s.	
LS-03M3-C185	HAISIL HL C18 5µm, 3cm x 0.3mm Microbore Column s.s.	
LS-03M5-C183	HAISIL HL C18 3µm, 3cm x 0.5mm Microbore Column s.s.	
LS-03M5-C185	HAISIL HL C18 5µm, 3cm x 0.5mm Microbore Column s.s.	
LS-0501-C183	HAISIL HL C18 3µm, 5cm x 1.0mm Microbore Column s.s.	
LS-0501-C185	HAISIL HL C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
LS-0510-C185	HAISIL HL C18 5µm, 5cm x 10mm Semi/Prep Column s.s.	
LS-0520-C185	HAISIL HL C18 5µm, 5cm x 20mm Semi/Prep Column s.s.	
LS-0521-C183	HAISIL HL C18 3µm, 5cm x 2.1mm Analytical Column s.s.	
LS-0521-C185	HAISIL HL C18 5µm, 5cm x 2.1mm Analytical Column s.s.	
LS-0530-C183	HAISIL HL C18 3µm, 5cm x 3.0mm Analytical Column s.s.	
LS-0530-C185	HAISIL HL C18 5µm, 5cm x 3.0mm Analytical Column s.s.	
LS-0546-C183	HAISIL HL C18 3µm, 5cm x 4.6mm Analytical Column s.s.	
LS-0546-C185	HAISIL HL C18 5µm, 5cm x 4.6mm Analytical Column s.s.	
LS-05M3-C183	HAISIL HL C18 3µm, 5cm x 0.3mm Microbore Column s.s.	
LS-05M3-C185	HAISIL HL C18 5µm, 5cm x 0.3mm Microbore Column s.s.	
LS-05M5-C183	HAISIL HL C18 3µm, 5cm x 0.5mm Microbore Column s.s.	
LS-05M5-C185	HAISIL HL C18 5µm, 5cm x 0.5mm Microbore Column s.s.	
LS-1001-C183	HAISIL HL C18 3µm, 10cm x 1.0mm Microbore Column s.s.	
LS-1001-C185	HAISIL HL C18 5µm, 10cm x 1.0mm Microbore Column s.s.	
LS-1010-C185	HAISIL HL C18 5µm, 10cm x 10mm Semi/Prep Column s.s.	
LS-1020-C185	HAISIL HL C18 5µm, 10cm x 20mm Semi/Prep Column s.s.	
LS-1021-C183	HAISIL HL C18 3µm, 10cm x 2.1mm Analytical Column s.s.	
LS-1021-C185	HAISIL HL C18 5µm, 10cm x 2.1mm Analytical Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
LS-1030-C183	HAISIL HL C18 3µm, 10cm x 3.0mm Analytical Column s.s.	
LS-1030-C185	HAISIL HL C18 5µm, 10cm x 3.0mm Analytical Column s.s.	
LS-1046-C183	HAISIL HL C18 3µm, 10cm x 4.6mm Analytical Column s.s.	
LS-1046-C185	HAISIL HL C18 5µm, 10cm x 4.6mm Analytical Column s.s.	
LS-10M3-C183	HAISIL HL C18 3µm, 10cm x 0.3mm Microbore Column s.s.	
LS-10M3-C185	HAISIL HL C18 5µm, 10cm x 0.3mm Microbore Column s.s.	
LS-10M5-C183	HAISIL HL C18 3µm, 10cm x 0.5mm Microbore Column s.s.	
LS-10M5-C185	HAISIL HL C18 5µm, 10cm x 0.5mm Microbore Column s.s.	
LS-1501-C183	HAISIL HL C18 3µm, 15cm x 1.0mm Microbore Column s.s.	
LS-1501-C185	HAISIL HL C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
LS-1510-C185	HAISIL HL C18 5µm, 15cm x 10mm Semi/Prep Column s.s.	
LS-1520-C185	HAISIL HL C18 5µm, 15cm x 20mm Semi/Prep Column s.s.	
LS-1521-C183	HAISIL HL C18 3µm, 15cm x 2.1mm Analytical Column s.s.	
LS-1521-C185	HAISIL HL C18 5µm, 15cm x 2.1mm Analytical Column s.s.	
LS-1530-C183	HAISIL HL C18 3µm, 15cm x 3.0mm Analytical Column s.s.	
LS-1530-C185	HAISIL HL C18 5µm, 15cm x 3.0mm Analytical Column s.s.	
LS-1546-C183	HAISIL HL C18 3µm, 15cm x 4.6mm Analytical Column s.s.	
LS-1546-C185	HAISIL HL C18 5µm, 15cm x 4.6mm Analytical Column s.s.	
LS-15M3-C185	HAISIL HL C18 5µm, 15cm x 0.3mm Microbore Column s.s.	
LS-15M5-C183	HAISIL HL C18 3µm, 15cm x 0.5mm Microbore Column s.s.	
LS-15M5-C185	HAISIL HL C18 5µm, 15cm x 0.5mm Microbore Column s.s.	
LS-2501-C185	HAISIL HL C18 5µm, 25cm x 1.0mm Microbore Column s.s.	
LS-2510-C185	HAISIL HL C18 5µm, 25cm x 10mm Semi/Prep Column s.s.	
LS-2520-C185	HAISIL HL C18 5µm, 25cm x 20mm Semi/Prep Column s.s.	
LS-2521-C183	HAISIL HL C18 3µm, 25cm x 2.1mm Analytical Column s.s.	
LS-2521-C185	HAISIL HL C18 5µm, 25cm x 2.1mm Analytical Column s.s.	
LS-2530-C183	HAISIL HL C18 3µm, 25cm x 3.0mm Analytical Column s.s.	
LS-2530-C185	HAISIL HL C18 5µm, 25cm x 3.0mm Analytical Column s.s.	
LS-2546-C183	HAISIL HL C18 3µm, 25cm x 4.6mm Analytical Column s.s.	
LS-2546-C185	HAISIL HL C18 5µm, 25cm x 4.6mm Analytical Column s.s.	
PC-0307-C185	PHALANX Capellini C18 5µm, 30 x 0.075mm Column	
PC-0315-C185	PHALANX Capellini C18 5µm, 30 x 0.150mm Column	
PC-0507-C185	PHALANX Capellini C18 5µm, 50 x 0.075mm Column	
PC-0515-C185	PHALANX Capellini C18 5µm, 50 x 0.150mm Column	
PC-1007-C185	PHALANX Capellini C18 5µm, 100 x 0.075mm Column	
PC-1015-C185	PHALANX Capellini C18 5µm, 100 x 0.150mm Column	
PC-1507-C185	PHALANX Capellini C18 5µm, 150 x 0.075mm Column	
PC-1515-C185	PHALANX Capellini C18 5µm, 150 x 0.150mm Column	
PC-2507-C185	PHALANX Capellini C18 5µm, 250 x 0.075mm Column	
PC-2515-C185	PHALANX Capellini C18 5µm, 250 x 0.150mm Column	
PF-0232-C185	PHALANX HAICART C18 5µm, 20 x 2.3mm 5/pk	
PF-0310-C185	PHALANX HAISpeed10 C18 5µm, 30 x 10mm	
PF-0432-C185	PHALANX HAICART C18 5µm, 40 x 3.2mm Cart. Column	
PF-0446-C185	PHALANX HAICART C18 5µm, 40 x 4.6mm Cart. Column	
PF-0510-C185	PHALANX HAISpeed10 C18 5µm, 50 x 10mm	
PF-0532-C185	PHALANX HAICART C18 5µm, 50 x 3.2mm Cart. Column	
PF-0546-C185	PHALANX HAICART C18 5µm, 50 x 4.6mm Cart. Column	
PF-1010-C185	PHALANX HAISpeed10 C18 5µm, 100 x 10mm	
PF-1032-C185	PHALANX HAICART C18 5µm, 100 x 3.2mm Cart. Column	
PF-1046-C185	PHALANX HAICART C18 5µm, 100 x 4.6mm Cart. Column	
PK-0221-C185	PHALANX HAIPEEK C18 5µm, 20 x 2.1mm	
PK-0421-C185	PHALANX HAIPEEK C18 5µm, 40 x 2.1mm	
PK-0521-C185	PHALANX HAIPEEK C18 5µm, 50 x 2.1mm	
PLC-0315-AXK5	POLYMERIC Capellini AXK 5µm, 30 x 0.150mm Column	
PLC-0315-CXK5	POLYMERIC Capellini RCXK 5µm, 30 x .150mm Column	
PLC-0315-RP13	POLYMERIC Capellini RPI 3µm, 30 x 0.150mm Column	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
PLC-0315-RP33	POLYMERIC Capellini RP3 3µm, 30 × 0.150mm Column	
PLC-0515-AXK5	POLYMERIC Capellini AXK 5µm, 50 × 0.150mm Column	
PLC-0515-CXK5	POLYMERIC Capellini CXK 5µm, 50 × 0.150mm Column	
PLC-0515-RP13	POLYMERIC Capellini RPI 3µm, 50 × 0.150mm Column	
PLC-0515-RP33	POLYMERIC Capellini RP3 3µm, 50 × 0.150mm Column	
PLC-1015-RP13	POLYMERIC Capellini RPI 3µm, 100 × 0.150mm Column	
PLC-1015-RP33	POLYMERIC Capellini RP3 3µm, 100 × 0.150mm Column	
PLC-1515-RP13	POLYMERIC Capellini RPI 3µm, 150 × 0.150mm Column	
PLC-1515-RP33	POLYMERIC Capellini RP3 3µm, 150 × 0.150mm Column	
PLC-2515-RP13	POLYMERIC Capellini RPI 3µm, 250 × 0.150mm Column	
PLC-2515-RP33	POLYMERIC Capellini RP3 3µm, 250 × 0.150mm Column	
PLF-0232-AXK5	POLYMERIC HAICART AXK 5µm, 20 × 3.2mm Cart. Column	
PLF-0232-CXK5	POLYMERIC HAICART CXK 5µm, 20 × 3.2mm Cart. Column	
PLF-0232-RP13	POLYMERIC HAICART RPI 3µm, 20 × 3.2mm Cart. Column	
PLF-0232-RP33	POLYMERIC HAICART RP3 3µm, 20 × 3.2mm Cart. Column	
PLK-0221-AXK5	POLYMERIC HAIPEEK AXK 5µm, 20 × 2.1mm Cart.Column	
PLK-0221-CXK5	POLYMERIC HAIPEEK CXK 5µm, 20 × 2.1mm Cart.Column	
PLK-0221-RP13	POLYMERIC HAIPEEK RPI 3µm, 20 × 2.1mm Cart.Column	
PLK-0221-RP33	POLYMERIC HAIPEEK RP3 3µm, 20 × 2.1mm Cart.Column	
PLK-0421-AXK5	POLYMERIC HAIPEEK AXK 5µm, 40 × 2.1mm Cart.Column	
PLK-0421-CXK5	POLYMERIC HAIPEEK CXK 5µm, 40 × 2.1mm Cart.Column	
PLK-0421-RP13	POLYMERIC HAIPEEK RPI 3µm, 40 × 2.1mm Cart.Column	
PLK-0421-RP33	POLYMERIC HAIPEEK RP3 3µm, 40 × 2.1mm Cart.Column	
PLK-0521-AXK5	POLYMERIC HAIPEEK AXK 5µm, 50 × 2.1mm Cart.Column	
PLK-0521-CXK5	POLYMERIC HAIPEEK CXK 5µm, 50 × 2.1mm Cart.Column	
PLK-0521-RP13	POLYMERIC HAIPEEK RPI 3µm, 50 × 2.1mm Cart.Column	
PLK-0521-RP33	POLYMERIC HAIPEEK RP3 3µm, 50 × 2.1mm Cart.Column	
PLP-0101-AXK5	POLYMERIC Piccolo AXK 5µm, 10 × 1.0mm Column	
PLP-0101-CXK5	POLYMERIC Piccolo CXK 5µm, 10 × 1.0mm Column	
PLP-0101-RP13	POLYMERIC Piccolo RPI 3µm, 10 × 1.0mm Column	
PLP-0101-RP33	POLYMERIC Piccolo RP3 3µm, 10 × 1.0mm Column	
PLP-M2M5-AXK5	POLYMERIC Piccolo AXK 5µm, 2.5 × 0.5mm Column	
PLP-M2M5-CXK5	POLYMERIC Piccolo CXK 5µm, 2.5 × 0.5mm Column	
PLP-M2M5-RP13	POLYMERIC Piccolo RPI 3µm, 2.5 × 0.5mm Column	
PLP-M2M5-RP33	POLYMERIC Piccolo RP3 3µm, 2.5 × 0.5mm Column	
PLP-M501-AXK5	POLYMERIC Piccolo AXK 5µm, 5 × 1.0mm Column	
PLP-M501-CXK5	POLYMERIC Piccolo CXK 5µm, 5 × 1.0mm Column	
PLP-M501-RP13	POLYMERIC Piccolo RPI 3µm, 5 × 1.0mm Column	
PLP-M501-RP33	POLYMERIC Piccolo RP3 3µm, 5 × 1.0mm Column	
PLR-0121-AXK5	POLYMERIC Sprite AXK 5µm, 10 × 2.1mm Column	
PLR-0121-CXK5	POLYMERIC Sprite CXK 5µm, 10 × 2.1mm Column	
PLR-0121-RP13	POLYMERIC Sprite RPI 3µm, 10 × 2.1mm Column	
PLR-0121-RP33	POLYMERIC Sprite RP3 3µm, 10 × 2.1mm Column	
PLR-0221-AXK5	POLYMERIC Sprite AXK 5µm, 20 × 2.1mm Column	
PLR-0221-CXK5	POLYMERIC Sprite CXK 5µm, 20 × 2.1mm Column	
PLR-0221-RP13	POLYMERIC Sprite RPI 3µm, 20 × 2.1mm Column	
PLR-0221-RP33	POLYMERIC Sprite RP3 3µm, 20 × 2.1mm Column	
PLR-0421-AXK5	POLYMERIC Sprite AXK 5µm, 40 × 2.1mm Column	
PLR-0421-CXK5	POLYMERIC Sprite CXK 5µm, 40 × 2.1mm Column	
PLR-0421-RP13	POLYMERIC Sprite RPI 3µm, 40 × 2.1mm Column	
PLR-0421-RP33	POLYMERIC Sprite RP3 3µm, 40 × 2.1mm Column	
PLS-0301-AXK5	POLYMERIC Microbore AXK 5µm, 30 × 1.0mm Column	
PLS-0301-CXK5	POLYMERIC Microbore CXK 5µm, 30 × 1.0mm Column	
PLS-0301-RP13	POLYMERIC Microbore RPI 3µm, 30 × 1.0mm Column	
PLS-0301-RP33	POLYMERIC Microbore RP3 3µm, 30 × 1.0mm Column	
PLS-03M3-AXK5	POLYMERIC Microbore AXK 5µm, 30 × 0.3mm Column	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
PLS-03M3-CXK5	POLYMERIC Microbore CXK 5µm, 30 X 0.3mm Column	
PLS-03M3-RP13	POLYMERIC Microbore RP1 3µm, 30 X 0.3mm Column	
PLS-03M3-RP33	POLYMERIC Microbore RP3 3µm, 30 X 0.3mm Column	
PLS-03M5-AXK5	POLYMERIC Microbore AXK 5µm, 30 x 0.5mm Column	
PLS-03M5-CXK5	POLYMERIC Microbore CXK 5µm, 30 x 0.5mm Column	
PLS-03M5-RP13	POLYMERIC Microbore RP1 3µm, 30 x 0.5mm Column	
PLS-03M5-RP33	POLYMERIC Microbore RP3 3µm, 30 x 0.5mm Column	
PLS-0501-AXK5	POLYMERIC Microbore AXK 5µm, 50 x 1.0mm Column	
PLS-0501-CXK5	POLYMERIC Microbore CXK 5µm, 50 x 1.0mm Column	
PLS-0501-RP13	POLYMERIC Microbore RP1 3µm, 50 x 1.0mm Column	
PLS-0501-RP33	POLYMERIC Microbore RP3 3µm, 50 x 1.0mm Column	
PLS-05M3-AXK5	POLYMERIC Microbore AXK 5µm, 50 x 0.3mm Column	
PLS-05M3-CXK5	POLYMERIC Microbore CXK 5µm, 50 x 0.3mm Column	
PLS-05M3-RP13	POLYMERIC Microbore RP1 3µm, 50 x 0.3mm Column	
PLS-05M3-RP33	POLYMERIC Microbore RP3 3µm, 50 x 0.3mm Column	
PLS-05M5-AXK5	POLYMERIC Microbore AXK 5µm, 50 x 0.5mm Column	
PLS-05M5-CXK5	POLYMERIC Microbore CXK 5µm, 50 x 0.5mm Column	
PLS-05M5-RP13	POLYMERIC Microbore RP1 3µm, 50 x 0.5mm Column	
PLS-05M5-RP33	POLYMERIC Microbore RP3 3µm, 50 x 0.5mm Column	
PLS-1001-AXK5	POLYMERIC Microbore AXK 5µm, 100 x 1.0mm Column	
PLS-1001-CXK5	POLYMERIC Microbore CXK 5µm, 100 x 1.0mm Column	
PLS-1001-RP13	POLYMERIC Microbore RP1 3µm, 100 x 1.0mm Column	
PLS-1001-RP33	POLYMERIC Microbore RP3 3µm, 100 x 1.0mm Column	
PLS-10M3-RP13	POLYMERIC Microbore RP1 3µm, 100 x 0.3MM Column	
PLS-10M3-RP33	POLYMERIC Microbore RP3 3µm, 100 x 0.3mm Column	
PLS-10M5-RP13	POLYMERIC Microbore RP1 3µm, 100 x 0.5mm Column	
PLS-10M5-RP33	POLYMERIC Microbore RP3 3µm, 100 x 0.5mm Column	
PLS-1501-AXK5	POLYMERIC Microbore AXK 5µm, 150 x 1.0mm Column	
PLS-1501-CXK5	POLYMERIC Microbore CXK 5µm, 150 x 1.0mm Column	
PLS-1501-RP13	POLYMERIC Microbore RP1 3µm, 100 x 1.0mm Column	
PLS-1501-RP33	POLYMERIC Microbore RP3 3µm, 150 x 1.0mm Column	
PLS-15M3-RP13	POLYMERIC Microbore RP1 3µm, 150 x 0.3mm Column	
PLS-15M3-RP33	POLYMERIC Microbore RP3 3µm, 150 x 0.3mm Column	
PLS-15M5-RP13	POLYMERIC Microbore RP1 3µm, 150 x 0.5mm Column	
PLS-15M5-RP33	POLYMERIC Microbore RP3 3µm, 150 x 0.5mm Column	
PLS-2501-RP13	POLYMERIC Microbore RP1 3µm, 250 x 1.0mm Column	
PLS-2501-RP33	POLYMERIC Microbore RP3 3µm, 250 x 1.0mm Column	
PLS-25M3-RP13	POLYMERIC Microbore RP1 3µm, 250 x 0.3mm Column	
PLS-25M3-RP33	POLYMERIC Microbore RP3 3µm, 250 x 0.3mm Column	
PLS-25M5-RP13	POLYMERIC Microbore RP1 3µm, 250 x 0.5mm Column	
PLS-25M5-RP33	POLYMERIC Microbore RP3 3µm, 250 x 0.5mm Column	
PP-0101-C185	PHALANX Piccolo C18 5µm, 10 x 1.0mm	
PP-M2M5-C185	PHALANX Piccolo C18 5µm, 2.5 x 0.5mm	
PP-M501-C185	PHALANX Piccolo C18 5µm, 5 x 1.0mm	
PR-0121-C185	PHALANX Sprite C18 5µm, 10 x 2.1mm	
PR-0221-C185	PHALANX Sprite C18 5µm, 20 x 2.1mm	
PR-0421-C185	PHALANX Sprite C18 5µm, 40 x 2.1mm	
PS-0301-C185	PHALANX C18 5µm, 30 x 0.1mm Column	
PS-0321-C185	PHALANX C18 5µm, 30 x 2.1mm Column s.s.	
PS-0330-C185	PHALANX C18 5µm, 30 x 3.0mm Column s.s.	
PS-0346-C185	PHALANX C18 5µm, 30 x 4.6mm Column s.s.	
PS-03M3-C185	PHALANX Microbore C18 5µm, 30 x 0.3mm	
PS-03M5-C185	PHALANX Microbore C18 5µm, 30 x 0.5mm	
PS-0501-C185	PHALANX Microbore C18 5µm, 50 x 1.0mm Column	
PS-0510-C185	PHALANX Semi-Prep C18 5µm, 50 x 10mm Column	
PS-0520-C185	PHALANX Prep C18 5µm, 50 x 20mm Column	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
PS-0521-C185	PHALANX C18 5µm, 50 × 2.1mm Column s.s.	
PS-0530-C185	PHALANX C18 5µm, 50 × 3.0mm Column s.s.	
PS-0546-C185	PHALANX C18 5µm, 50 × 4.6mm Column s.s.	
PS-05M3-C185	PHALANX Microbore 5µm, 5µm 50 × 0.3mm Column	
PS-05M5-C185	PHALANX Microbore 5µm, 5µm 50 × 0.5mm Column	
PS-1001-C185	PHALANX Microbore C18 5µm, 100 × 1.0mm Column	
PS-1010-C185	PHALANX Semi-Prep C18 5µm, 100 × 10mm Column	
PS-1020-C185	PHALANX Prep C18 5µm, 100 × 20mm Column	
PS-1021-C185	PHALANX C18 5µm, 100 × 2.1mm Column s.s.	
PS-1030-C185	PHALANX C18 5µm, 100 × 3.0mm Column s.s.	
PS-1046-C185	PHALANX C18 5µm, 100 × 4.6mm Column s.s.	
PS-10M3-C185	PHALANX Microbore C18 5µm, 100 × 0.3mm Column	
PS-10M5-C185	PHALANX Microbore C18 5µm, 100 × 0.5mm Column	
PS-1501-C185	PHALANX Microbore C18 5µm, 150 × 1.0mm Column	
PS-1510-C185	PHALANX Semi-Prep C18 5µm, 150 × 10mm Column	
PS-1520-C185	PHALANX Prep C18 5µm, 150 × 20mm Column	
PS-1521-C185	PHALANX C18 5µm, 150 × 2.1mm Column s.s.	
PS-1530-C185	PHALANX C18 5µm, 150 × 3.0mm Column s.s.	
PS-1546-C185	PHALANX C18 5µm, 150 × 4.6mm Column s.s.	
PS-15M3-C185	PHALANX Microbore 5µm, 5µm 150 × 0.3mm Column	
PS-15M5-C185	PHALANX Microbore C18 5µm, 150 × 0.5mm Column	
PS-2501-C185	PHALANX Microbore C18 5µm, 250 × 1.0mm Column	
PS-2510-C185	PHALANX Semi-Prep C18 5µm, 250 × 10mm Column	
PS-2520-C185	PHALANX Prep C18 5µm, 250 × 20mm Column	
PS-2521-C185	PHALANX C18 5µm, 250 × 2.1mm Column s.s.	
PS-2530-C185	PHALANX C18 5µm, 250 × 3.0mm Column s.s.	
PS-2546-C185	PHALANX C18 5µm, 250 × 4.6mm Column s.s.	
PS-25M3-C185	PHALANX Microbore C18 5µm, 250 × 0.3mm Column	
PS-25M5-C185	PHALANX Microbore C18 5µm, 250 × 0.5mm Column	
RP-0101-C045	PROTEO300 Piccolo C4 5µm, 10 × 1.0mm	
RP-0101-C185	PROTEO300 Piccolo C18 5µm, 10 × 1.0mm	
RP-M2M5-C045	PROTEO300 Piccolo C4 5µm, 2.5 × 0.5mm	
RP-M2M5-C185	PROTEO300 Piccolo C18 5µm, 2.5 × 0.5mm	
RP-M501-C045	PROTEO300 Piccolo C4 5µm, 5 × 1.0mm	
RP-M501-C185	PROTEO300 Piccolo C18 5µm, 5 × 1.0mm	
RR-0121-C045	PROTEO300 Sprite C4 5µm, 10 × 2.1mm	
RR-0121-C185	PROTEO300 Sprite C18 5µm, 10 × 2.1mm	
RR-0221-C045	PROTEO300 Sprite C4 5µm, 20 × 2.1mm	
RR-0221-C185	PROTEO300 Sprite C18 5µm, 20 × 2.1mm	
RR-0421-C045	PROTEO300 Sprite C4 5µm, 40 × 2.1mm	
RR-0421-C185	PROTEO300 Sprite C18 5µm, 40 × 2.1mm	
RS-0301-C045	PROTEO300 C4 5µm, 30 × 0.1mm Column	
RS-0301-C185	PROTEO300 C18 5µm, 30 × 0.1mm Column	
RS-0321-C045	PROTEO300 C4 5µm, 30 × 2.1mm Column s.s.	
RS-0321-C185	PROTEO300 C18 5µm, 30 × 2.1mm Column s.s.	
RS-0330-C045	PROTEO300 C4 5µm, 30 × 3.0mm Column s.s.	
RS-0330-C185	PROTEO300 C18 5µm, 30 × 3.0mm Column s.s.	
RS-0346-C045	PROTEO300 C4 5µm, 30 × 4.6mm Column s.s.	
RS-0346-C185	PROTEO300 C18 5µm, 30 × 4.6mm Column s.s.	
RS-03M3-C045	PROTEO300 Microbore C4 5µm, 30 × 0.3mm	
RS-03M3-C185	PROTEO300 Microbore C18 5µm, 30 × 0.3mm	
RS-03M5-C045	PROTEO300 Microbore C4 5µm, 30 × 0.5mm	
RS-03M5-C185	PROTEO300 Microbore C18 5µm, 30 × 0.5mm	
RS-0501-C045	PROTEO300 Microbore C4 5µm, 50 × 1.0mm Column	
RS-0501-C185	PROTEO300 Microbore C18 5µm, 50 × 1.0mm Column	
RS-0510-C041	PROTEO300 Semi-Prep C4 5µm, 50 × 10mm Column	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
RS-0510-C045	PROTEO300 Semi-Prep C4 5µm, 50 × 10mm Column	
RS-0510-C181	PROTEO300 Semi-Prep C18 5µm, 50 × 10mm Column	
RS-0510-C185	PROTEO300 Semi-Prep C18 5µm, 50 × 10mm Column	
RS-0520-C041	PROTEO300 Prep C4 10µm, 50 × 20mm Column	
RS-0520-C045	PROTEO300 Prep C4 5µm, 50 × 20mm Column	
RS-0520-C181	PROTEO300 Prep C18 10µm, 50 × 20mm Column	
RS-0520-C185	PROTEO300 Prep C18 5µm, 50 × 20mm Column	
RS-0521-C045	PROTEO300 C4 5µm, 50 × 2.1mm Column s.s.	
RS-0521-C185	PROTEO300 C18 5µm, 50 × 2.1mm Column s.s.	
RS-0530-C045	PROTEO300 C4 5µm, 50 × 3.0mm Column s.s.	
RS-0530-C185	PROTEO300 C18 5µm, 50 × 3.0mm Column s.s.	
RS-0546-C045	PROTEO300 C4 5µm, 50 × 4.6mm Column s.s.	
RS-0546-C185	PROTEO300 C18 5µm, 50 × 4.6mm Column s.s.	
RS-05M3-C045	PROTEO300 Microbore 5µm, 5µm 50 × 0.3mm Column	
RS-05M3-C185	PROTEO300 Microbore 5µm, 5µm 50 × 0.3mm Column	
RS-05M5-C045	PROTEO300 Microbore 5µm, 5µm 50 × 0.5mm Column	
RS-05M5-C185	PROTEO300 Microbore 5µm, 5µm 50 × 0.5mm Column	
RS-1001-C045	PROTEO300 Microbore C4 5µm, 100 × 1.0mm Column	
RS-1001-C185	PROTEO300 Microbore C18 5µm, 100 × 1.0mm Column	
RS-1010-C041	PROTEO300 Semi-Prep C4 10µm, 100 × 10mm Column	
RS-1010-C045	PROTEO300 Semi-Prep C4 5µm, 100 × 10mm Column	
RS-1010-C181	PROTEO300 Semi-Prep C18 10µm, 100 × 10mm Column	
RS-1010-C185	PROTEO300 Semi-Prep C18 5µm, 100 × 10mm Column	
RS-1020-C041	PROTEO300 Prep C4 10µm, 100 × 20mm Column	
RS-1020-C045	PROTEO300 Prep C4 5µm, 100 × 20mm Column	
RS-1020-C181	PROTEO300 Prep C18 10µm, 100 × 20mm Column	
RS-1020-C185	PROTEO300 Prep C18 5µm, 100 × 20mm Column	
RS-1021-C045	PROTEO300 C4 5µm, 100 × 2.1mm Column s.s.	
RS-1021-C185	PROTEO300 C18 5µm, 100 × 2.1mm Column s.s.	
RS-1030-C045	PROTEO300 C4 5µm, 100 × 3.0mm Column s.s.	
RS-1030-C185	PROTEO300 C18 5µm, 100 × 3.0mm Column s.s.	
RS-1046-C045	PROTEO300 C4 5µm, 100 × 4.6mm Column s.s.	
RS-1046-C185	PROTEO300 C18 5µm, 100 × 4.6mm Column s.s.	
RS-10M3-C045	PROTEO300 Microbore C4 5µm, 100 × 0.3mm Column	
RS-10M3-C185	PROTEO300 Microbore C18 5µm, 100 × 0.3mm Column	
RS-10M5-C045	PROTEO300 Microbore C4 5µm, 100 × 0.5mm Column	
RS-10M5-C185	PROTEO300 Microbore C18 5µm, 100 × 0.5mm Column	
RS-1501-C045	PROTEO300 Microbore C4 5µm, 150 × 1.0mm Column	
RS-1501-C185	PROTEO300 Microbore C18 5µm, 150 × 1.0mm Column	
RS-1510-C041	PROTEO300 Semi-Prep C4 10µm, 150 × 10mm Column	
RS-1510-C045	PROTEO300 Semi-Prep C4 5µm, 150 × 10mm Column	
RS-1510-C181	PROTEO300 Semi-Prep C18 10µm, 150 × 10mm Column	
RS-1510-C185	PROTEO300 Semi-Prep C18 5µm, 150 × 10mm Column	
RS-1520-C041	PROTEO300 Prep C4 10µm, 150 × 20mm Column	
RS-1520-C045	PROTEO300 Prep C4 5µm, 150 × 20mm Column	
RS-1520-C181	PROTEO300 Prep C18 10µm, 150 × 20mm Column	
RS-1520-C185	PROTEO300 Prep C18 5µm, 150 × 20mm Column	
RS-1521-C045	PROTEO300 C4 5µm, 150 × 2.1mm Column s.s.	
RS-1521-C185	PROTEO300 C18 5µm, 150 × 2.1mm Column s.s.	
RS-1530-C045	PROTEO300 C4 5µm, 150 × 3.0mm Column s.s.	
RS-1530-C185	PROTEO300 C18 5µm, 150 × 3.0mm Column s.s.	
RS-1546-C045	PROTEO300 C4 5µm, 150 × 4.6mm Column s.s.	
RS-1546-C185	PROTEO300 C18 5µm, 150 × 4.6mm Column s.s.	
RS-15M3-C045	PROTEO300 Microbore 5µm, 5µm 150 × 0.3mm Column	
RS-15M3-C185	PROTEO300 Microbore 5µm, 5µm 150 × 0.3mm Column	
RS-15M5-C045	PROTEO300 Microbore C4 5µm, 150 × 0.5mm Column	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
RS-15M5-C185	PROTEO300 Microbore C18 5µm, 150 x 0.5mm Column	
RS-2501-C045	PROTEO300 Microbore C4 5µm, 250 x 1.0mm Column	
RS-2501-C185	PROTEO300 Microbore C18 5µm, 250 x 1.0mm Column	
RS-2510-C041	PROTEO300 Semi-Prep C4 10µm, 250 x 10mm Column	
RS-2510-C045	PROTEO300 Semi-Prep C4 5µm, 250 x 10mm Column	
RS-2510-C181	PROTEO300 Semi-Prep C18 10µm, 250 x 10mm Column	
RS-2510-C185	PROTEO300 Semi-Prep C18 5µm, 250 x 10mm Column	
RS-2520-C041	PROTEO300 Prep C4 10µm, 250 x 20mm Column	
RS-2520-C045	PROTEO300 Prep C4 5µm, 250 x 20mm Column	
RS-2520-C181	PROTEO300 Prep C18 10µm, 250 x 20mm Column	
RS-2520-C185	PROTEO300 Prep C18 5µm, 250 x 20mm Column	
RS-2521-C045	PROTEO300 C4 5µm, 250 x 2.1mm Column s.s.	
RS-2521-C185	PROTEO300 C18 5µm, 250 x 2.1mm Column s.s.	
RS-2530-C045	PROTEO300 C4 5µm, 250 x 3.0mm Column s.s.	
RS-2530-C185	PROTEO300 C18 5µm, 250 x 3.0mm Column s.s.	
RS-2546-C045	PROTEO300 C4 5µm, 250 x 4.6mm Column s.s.	
RS-2546-C185	PROTEO300 C18 5µm, 250 x 4.6mm Column s.s.	
RS-25M3-C045	PROTEO300 Microbore C4 5µm, 250 x 0.3mm Column	
RS-25M3-C185	PROTEO300 Microbore C18 5µm, 250 x 0.3mm Column	
RS-25M5-C045	PROTEO300 Microbore C4 5µm, 250 x 0.5mm Column	
RS-25M5-C185	PROTEO300 Microbore C18 5µm, 250 x 0.5mm Column	
T108	TARGA C8 5µm, Guard Cartridge 5/pkg	
T118	TARGA C18 5µm, Guard Cartridge s.s. 5/pkg	
TC-0507-C185	TARGA Capellini C18 5µm, 50 x 0.150mm column	
TC-0515-C185	TARGA Capellini C18 5µm, 50 x 0.150mm column	
TC-1007-C185	TARGA Capellini C18 5µm, 100 x 0.075mm column	
TC-1015-C185	TARGA Capellini C18 5µm, 100 x 0.150mm column	
TC-1507-C185	TARGA Capellini C18 5µm, 150 x 0.075mm column	
TC-1515-C185	TARGA Capellini C18 5µm, 150 x 0.150mm column	
TC-2507-C185	TARGA Capellini C18 5µm, 250 x 0.075mm column	
TC-2515-C185	TARGA Capellini C18 5µm, 250 x 0.150mm column	
TF-0220-C181	TARGA C18 10µm, 2.5cm x 20mm Cartridge s.s. 2/pkg	
TF-0310-C085	TARGA C8 5µm, 3cm x 10mm Cartridge Column s.s. 1/pkg	
TF-0310-C181	TARGA C18 10µm, 3cm x 10m Cartridge Column s.s.	
TF-0310-C185	TARGA C18 5µm, 3cm x 10m Cartridge Column s.s.	
TF-0432-C085	TARGA C8 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
TF-0432-C185	TARGA C18 5µm, 4cm x 3.2mm Cartridge s.s. 2/pkg	
TF-0446-C085	TARGA C8 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
TF-0446-C185	TARGA C18 5µm, 4cm x 4.6mm Cartridge s.s. 2/pkg	
TF-0510-C085	TARGA C8 5µm, 5cm x 10mm Cartridgr Column s.s.	
TF-0510-C181	TARGA C18 10µm, 5cm x 10mm Cartridge s.s. 2/pkg	
TF-0510-C185	TARGA C18 5µm, 5cm x 10mm Cartridge s.s.	
TF-0532-C085	TARGA C8 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
TF-0532-C185	TARGA C18 5µm, 5cm x 3.2mm Cartridge s.s. 2/pkg	
TF-0546-C085	TARGA C8 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
TF-0546-C185	TARGA C18 5µm, 5cm x 4.6mm Cartridge s.s. 2/pkg	
TF-1010-C085	TARGA C8 5µm, 10cm x 10mm Cartridge Column s.s.	
TF-1010-C181	TARGA C18 10µm, 10cm x 10mm Cartridge Column s.s.	
TF-1010-C185	TARGA C18 5µm, 10cm x 10mm Cartridge s.s.	
TF-1032-C085	TARGA C8 5µm, 10cm x 3.2mm Cartridge s.s.	
TF-1032-C185	TARGA C18 5µm, 10cm x 3.2mm Cartridge s.s.	
TF-1046-C085	TARGA C8 5µm, 10cm x 4.6mm Cartridge s.s.	
TF-1046-C185	TARGA C18 5µm, 10cm x 4.6mm Cartridge s.s.	
TK-0221-C085	TARGA C8 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	
TK-0221-C183	TARGA C18 3µm, 2cm x 2.1mm Guard Cartridge Peek	
TK-0221-C185	TARGA C18 5µm, 2cm x 2.1mm Cartridge PEEK 5/pkg	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
TK-0421-C085	TARGA C8 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
TK-0421-C183	TARGA C18 3µm, 4cm x 2.1mm Analytical Cartridge PEEK	
TK-0421-C185	TARGA C18 5µm, 4cm x 2.1mm Cartridge PEEK 2/pkg	
TK-0521-C085	TARGA C8 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
TK-0521-C183	TARGA C18 3µm, 5cm x 2.1mm Analytical Cartridge PEEK	
TK-0521-C185	TARGA C18 5µm, 5cm x 2.1mm Cartridge PEEK 2/pkg	
TP-0101-C185	TARGA Piccolo C18 5µm, 10 x 1.0mm column	
TP-M2M5-C185	TARGA Piccolo C18 5µm, 2.5 x 0.5mm column	
TP-M501-C185	TARGA Piccolo C18 5µm, 5 x 1.0mm column	
TR-0121-C185	TARGA Sprite C18 5µm, 10 x 2.1mm column	
TR-0221-C185	TARGA Sprite C18 5µm, 20 x 2.1mm column	
TR-0421-C185	TARGA Sprite C18 5µm, 40 x 2.1mm column	
TS-0301-C085	TARGA C8 5µm, 3cm x 1.0mm Microbore Column s.s.	
TS-0301-C183	TARGA C18 3µm, 3cm x 1.0mm Microbore Column s.s.	
TS-0301-C185	TARGA C18 5µm, 3cm x 1.0mm Microbore Column s.s.	
TS-0321-C085	TARGA C8 5µm, 3cm x 2.1mm Analytical Column s.s.	
TS-0321-C183	TARGA C18 3µm, 3cm x 2.1mm Analytical Column s.s.	
TS-0321-C185	TARGA C18 5µm, 3cm x 2.1mm Analytical Column s.s.	
TS-0330-C085	TARGA C8 5µm, 3cm x 3.0mm Analytical Column s.s.	
TS-0330-C183	TARGA C18 3µm, 3cm x 3.0mm Analytical Column s.s.	
TS-0330-C185	TARGA C18 5µm, 3cm x 3.0mm Analytical Column s.s.	
TS-0346-C085	TARGA C8 5µm, 3cm x 4.6mm Analytical Column s.s.	
TS-0346-C183	TARGA C18 3µm, 3cm x 4.6mm Analytical Column s.s.	
TS-0346-C185	TARGA C18 5µm, 3cm x 4.6mm Analytical Column s.s.	
TS-03M3-C085	TARGA C8 5µm, 3cm x 0.3mm Microbore Column s.s.	
TS-03M3-C183	TARGA C18 3µm, 3cm x 0.3mm Microbore Column s.s.	
TS-03M3-C185	TARGA C18 5µm, 3cm x 0.3mm Microbore Column s.s.	
TS-03M5-C085	TARGA C8 5µm, 3cm x 0.5mm Microbore Column s.s.	
TS-03M5-C183	TARGA C18 3µm, 3cm x 0.5mm Microbore Column s.s.	
TS-03M5-C185	TARGA C18 5µm, 3cm x 0.5mm Microbore Column s.s.	
TS-0501-C085	TARGA C8 5µm, 5cm x 1.0mm Microbore Column s.s.	
TS-0501-C183	TARGA C18 3µm, 5cm x 1.0mm Microbore Column s.s.	
TS-0501-C185	TARGA C18 5µm, 5cm x 1.0mm Microbore Column s.s.	
TS-0510-C085	TARGA C8 5µm, 5cm x 10mm Semi/Prep Column s.s.	
TS-0510-C181	TARGA C18 10µm, 5cm x 10mm Semi/Prep Column s.s.	
TS-0510-C185	TARGA C18 5µm, 5cm x 10mm Semi/Prep Column s.s.	
TS-0520-C085	TARGA C8 5µm, 5cm x 20mm Semi/Prep Column s.s.	
TS-0520-C181	TARGA C18 10µm, 5cm x 20mm Semi/Prep Column s.s.	
TS-0520-C185	TARGA C18 5µm, 5cm x 20mm Semi/Prep Column s.s.	
TS-0521-C085	TARGA C8 5µm, 5cm x 2.1mm Analytical Column s.s.	
TS-0521-C183	TARGA C18 3µm, 5cm x 2.1mm Analytical Column s.s.	
TS-0521-C185	TARGA C18 5µm, 5cm x 2.1mm Analytical Column s.s.	
TS-0530-C085	TARGA C8 5µm, 5cm x 3.0mm Analytical Column s.s.	
TS-0530-C183	TARGA C18 3µm, 5cm x 3.0mm Analytical Column s.s.	
TS-0530-C185	TARGA C18 5µm, 5cm x 3.0mm Analytical Column s.s.	
TS-0546-C085	TARGA C8 5µm, 5cm x 4.6mm Analytical Column s.s.	
TS-0546-C183	TARGA C18 3µm, 5cm x 4.6mm Analytical Column s.s.	
TS-0546-C185	TARGA C18 5µm, 5cm x 4.6mm Analytical Column s.s.	
TS-05M3-C085	TARGA C8 5µm, 5cm x 0.3mm Microbore Column s.s.	
TS-05M3-C183	TARGA C18 5µm, 5cm x 0.3mm Microbore Column s.s.	
TS-05M3-C185	TARGA C18 5µm, 5cm x 0.3mm Microbore Column s.s.	
TS-05M5-C085	TARGA C8 5µm, 5cm x 0.5mm Microbore Column s.s.	
TS-05M5-C183	TARGA C18 3µm, 5cm x 0.5mm Microbore Column s.s.	
TS-05M5-C185	TARGA C18 5µm, 5cm x 0.5mm Microbore Column s.s.	
TS-1001-C085	TARGA C8 5µm, 10cm x 1.0mm Microbore Column s.s.	
TS-1001-C183	TARGA C18 3µm, 10cm x 1.0mm Microbore Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
TS-1001-C185	TARGA C18 5µm, 10cm x 1.0mm Microbore Column s.s.	
TS-1010-C085	TARGA C8 5µm, 10cm x 10mm Semi/Prep Column s.s.	
TS-1010-C181	TARGA C18 10µm, 10cm x 10mm Semi/Prep Column s.s.	
TS-1010-C185	TARGA C18 5µm, 10cm x 10mm Semi/Prep Column s.s.	
TS-1020-C085	TARGA C8 5µm, 10cm x 10mm Semi/Prep Column s.s.	
TS-1020-C181	TARGA C18 10µm, 10cm x 20mm Semi/Prep Column s.s.	
TS-1020-C185	TARGA C18 5µm, 10cm x 20mm Semi/Prep Column s.s.	
TS-1021-C085	TARGA C8 5µm, 10cm x 2.1mm Analytical Column s.s.	
TS-1021-C183	TARGA C18 3µm, 10cm x 2.1mm Analytical Column s.s.	
TS-1021-C185	TARGA C18 5µm, 10cm x 2.1mm Analytical Column s.s.	
TS-1030-C085	TARGA C8 5µm, 10cm x 3.0mm Analytical Column s.s.	
TS-1030-C183	TARGA C18 3µm, 10cm x 3.0mm Analytical Column s.s.	
TS-1030-C185	TARGA C18 5µm, 10cm x 3.0mm Analytical Column s.s.	
TS-1046-C085	TARGA C8 5µm, 10cm x 4.6mm Analytical Column s.s.	
TS-1046-C183	TARGA C18 3µm, 10cm x 4.6mm Analytical Column s.s.	
TS-1046-C185	TARGA C18 5µm, 10cm x 4.6mm Analytical Column s.s.	
TS-10M3-C085	TARGA C8 5µm, 10cm x 0.3mm Microbore Column s.s.	
TS-10M3-C183	TARGA C18 3µm, 10cm x 0.3mm Microbore Column s.s.	
TS-10M3-C185	TARGA C18 5µm, 10cm x 0.3mm Microbore Column s.s.	
TS-10M5-C085	TARGA C8 5µm, 10cm x 0.5mm Microbore Column s.s.	
TS-10M5-C183	TARGA C18 3µm, 10cm x 0.5mm Microbore Column s.s.	
TS-10M5-C185	TARGA C18 5µm, 10cm x 0.5mm Microbore Column s.s.	
TS-1501-C085	TARGA C8 5µm, 15cm x 1.0mm Microbore Column s.s.	
TS-1501-C183	TARGA C18 5µm, 15cm x 1.0mm Microbore Column s.s.	
TS-1501-C185	TARGA C18 5µm, 15cm x 1.0mm Microbore Column s.s.	
TS-1510-C085	TARGA C8 5µm, 15cm x 10mm Semi/Prep Column s.s.	
TS-1510-C181	TARGA C18 10µm, 15cm x 10mm Semi/Prep Column s.s.	
TS-1510-C185	TARGA C18 5µm, 15cm x 10mm Semi/Prep Column s.s.	
TS-1520-C085	TARGA C8 5µm, 15cm x 10mm Semi/Prep Column s.s.	
TS-1520-C181	TARGA C18 10µm, 15cm x 20mm Semi/Prep Column s.s.	
TS-1520-C185	TARGA C18 5µm, 15cm x 20mm Semi/Prep Column s.s.	
TS-1521-C085	TARGA C8 5µm, 15cm x 2.1mm Analytical Column s.s.	
TS-1521-C183	TARGA C18 3µm, 15cm x 2.1mm Analytical Column s.s.	
TS-1521-C185	TARGA C18 5µm, 15cm x 2.1mm Analytical Column s.s.	
TS-1530-C085	TARGA C8 5µm, 15cm x 3.0mm Analytical Column s.s.	
TS-1530-C183	TARGA C18 3µm, 15cm x 3.0mm Analytical Column s.s.	
TS-1530-C185	TARGA C18 5µm, 15cm x 3.0mm Analytical Column s.s.	
TS-1546-C085	TARGA C8 5µm, 15cm x 4.6mm Analytical Column s.s.	
TS-1546-C183	TARGA C18 3µm, 15cm x 4.6mm Analytical Column s.s.	
TS-1546-C185	TARGA C18 5µm, 15cm x 4.6mm Analytical Column s.s.	
TS-15M5-C085	TARGA C8 5µm, 15cm x 0.5mm Microbore Column s.s.	
TS-15M5-C183	TARGA C18 3µm, 15cm x 0.5mm Microbore Column s.s.	
TS-15M5-C185	TARGA C18 5µm, 15cm x 0.5mm Microbore Column s.s.	
TS-2501-C085	TARGA C8 5µm, 25cm x 1.0mm Microbore Column s.s.	
TS-2501-C183	TARGA C18 3µm, 25cm x 1.0mm Microbore Column s.s.	
TS-2501-C185	TARGA C18 5µm, 25cm x 1.0mm Microbore Column s.s.	
TS-2510-C085	TARGA C8 5µm, 25cm x 10mm Semi/Prep Column s.s.	
TS-2510-C181	TARGA C18 10µm, 25cm x 10mm Semi/Prep Column s.s.	
TS-2510-C185	TARGA C18 5µm, 25cm x 10mm Semi/Prep Column s.s.	
TS-2520-C085	TARGA C8 5µm, 25cm x 10mm Semi/Prep Column s.s.	
TS-2520-C181	TARGA C18 10µm, 25cm x 20mm Semi/Prep Column s.s.	
TS-2520-C185	TARGA C18 5µm, 25cm x 20mm Semi/Prep Column s.s.	
TS-2521-C085	TARGA C8 5µm, 25cm x 2.1mm Analytical Column s.s.	
TS-2521-C183	TARGA C18 3µm, 25cm x 2.1mm Analytical Column s.s.	
TS-2521-C185	TARGA C18 5µm, 25cm x 2.1mm Analytical Column s.s.	
TS-2530-C085	TARGA C8 5µm, 25cm x 3.0mm Analytical Column s.s.	

HIGGINS ANALYTICAL 2006 FACTORY PRICE LIST

<i>Part Number</i>	<i>Description</i>	<i>Price</i>
TS-2530-C183	TARGA C18 3µm, 25cm x 3.0mm Analytical Column s.s.	
TS-2530-C185	TARGA C18 5µm, 25cm x 3.0mm Analytical Column s.s.	
TS-2546-C085	TARGA C8 5µm, 25cm x 4.6mm Analytical Column s.s.	
TS-2546-C183	TARGA C18 3µm, 25cm x 4.6mm Analytical Column s.s.	
TS-2546-C185	TARGA C18 5µm, 25cm x 4.6mm Analytical Column s.s.	
VP-0101-C185	HEAVY Piccolo C18 5µm, 10 x 1.0mm	
VP-M2M5-C185	HEAVY Piccolo C18 5µm, 2.5 x 0.5mm	
VP-M501-C185	HEAVY Piccolo C18 5µm, 5 x 1.0mm	
VR-0121-C185	HEAVY Sprite C18 5µm, 10 x 2.1mm	
VR-0221-C185	HEAVY Sprite C18 5µm, 20 x 2.1mm	
VR-0421-C185	HEAVY Sprite C18 5µm, 40 x 2.1mm	
VS-0301-C185	HEAVY C18 5µm, 30 x 0.1mm Column	
VS-0321-C185	HEAVY C18 5µm, 30 x 2.1mm Column s.s.	
VS-0330-C185	HEAVY C18 5µm, 30 x 3.0mm Column s.s.	
VS-0346-C185	HEAVY C18 5µm, 30 x 4.6mm Column s.s.	
VS-03M3-C185	HEAVY Microbore C18 5µm, 30 x 0.3mm	
VS-03M5-C185	HEAVY Microbore C18 5µm, 30 x 0.5mm	
VS-0501-C185	HEAVY Microbore C18 5µm, 50 x 1.0mm Column	
VS-0510-C185	HEAVY Semi-Prep C18 5µm, 50 x 10mm Column	
VS-0520-C185	HEAVY Prep C18 5µm, 50 x 20mm Column	
VS-0521-C185	HEAVY C18 5µm, 50 x 2.1mm Column s.s.	
VS-0530-C185	HEAVY C18 5µm, 50 x 3.0mm Column s.s.	
VS-0546-C185	HEAVY C18 5µm, 50 x 4.6mm Column s.s.	
VS-05M3-C185	HEAVY Microbore 5µm, 5µm 50 x 0.3mm Column	
VS-05M5-C185	HEAVY Microbore 5µm, 5µm 50 x 0.5mm Column	
VS-1001-C185	HEAVY Microbore C18 5µm, 100 x 1.0mm Column	
VS-1010-C185	HEAVY Semi-Prep C18 5µm, 100 x 10mm Column	
VS-1020-C185	HEAVY Prep C18 5µm, 100 x 20mm Column	
VS-1021-C185	HEAVY C18 5µm, 100 x 2.1mm Column s.s.	
VS-1030-C185	HEAVY C18 5µm, 100 x 3.0mm Column s.s.	
VS-1046-C185	HEAVY C18 5µm, 100 x 4.6mm Column s.s.	
VS-10M3-C185	HEAVY Microbore C18 5µm, 100 x 0.3mm Column	
VS-10M5-C185	HEAVY Microbore C18 5µm, 100 x 0.5mm Column	
VS-1501-C185	HEAVY Microbore C18 5µm, 150 x 1.0mm Column	
VS-1510-C185	HEAVY Semi-Prep C18 5µm, 150 x 10mm Column	
VS-1520-C185	HEAVY Prep C18 5µm, 150 x 20mm Column	
VS-1521-C185	HEAVY C18 5µm, 150 x 2.1mm Column s.s.	
VS-1530-C185	HEAVY C18 5µm, 150 x 3.0mm Column s.s.	
VS-1546-C185	HEAVY C18 5µm, 150 x 4.6mm Column s.s.	
VS-15M3-C185	HEAVY Microbore 5µm, 5µm 150 x 0.3mm Column	
VS-15M5-C185	HEAVY Microbore C18 5µm, 150 x 0.5mm Column	
VS-2501-C185	HEAVY Microbore C18 5µm, 250 x 1.0mm Column	
VS-2510-C185	HEAVY Semi-Prep C18 5µm, 250 x 10mm Column	
VS-2520-C185	HEAVY Prep C18 5µm, 250 x 20mm Column	
VS-2521-C185	HEAVY C18 5µm, 250 x 2.1mm Column s.s.	
VS-2530-C185	HEAVY C18 5µm, 250 x 3.0mm Column s.s.	
VS-2546-C185	HEAVY C18 5µm, 250 x 4.6mm Column s.s.	
VS-25M3-C185	HEAVY Microbore C18 5µm, 250 x 0.3mm Column	
VS-25M5-C185	HEAVY Microbore C18 5µm, 250 x 0.5mm Column	
WK-GUARD-FM	HAIPPEEK Guard Holder Kit (female/male) for Waters Columns	

HIGGINS ANALYTICAL PRODUCTS MAY BE PURCHASED AND ARE SOLD ONLY ON THE FOLLOWING WARRANTY TERMS AND CONDITIONS. ADDITIONAL OR DIFFERENT TERMS, WHETHER IN AN ORDER, ACKNOWLEDGMENT, INVOICE OR OTHER DOCUMENT, SHALL NOT BECOME PART OF OR VARY THE FOLLOWING TERMS UNLESS AGREED TO IN WRITING BY HIGGINS ANALYTICAL.

LIMITED WARRANTY

Higgins Analytical products are manufactured and quality controlled to the highest industry standards, and Higgins Analytical warrants that the products, when properly applied in accordance with Higgins Analytical's standards therefor, will perform for a period of thirty (30) days after delivery substantially as described in Higgins Analytical's then current specifications for the product. If a customer believes a product does not comply with such warranty, the customer shall, within the warranty period, report the purported defect to Higgins Analytical and cooperate fully with Higgins Analytical in determining the conditions and cause of any purported defect and, thereby, obtain a returned goods authorization from Higgins Analytical or the place of purchase. If Higgins Analytical determines the alleged defect occurred as a result of misuse, neglect, improper installation, repair, alteration, application or accident, rapid depressurization, removal of column end fittings, unusual physical or electrical stress or improper handling procedure, the warranty shall be void and ineffective.

EXPRESS LIMITATIONS

HIGGINS ANALYTICAL PRODUCTS ARE MANUFACTURED FOR USE IN RESEARCH IN CHEMICAL MANUFACTURING LABORATORIES ONLY AND ARE NOT INTENDED FOR USE IN CLINICAL OR DIAGNOSTIC APPLICATIONS. THE HAI SPEED 10 AND HAI SPEED 20 CARTRIDGE HOLDERS ARE RATED FOR USE UP TO 3,000 PSI (200 BAR) WITH FLUIDS THAT ARE LIQUID AT STP (STANDARD TEMPERATURE AND PRESSURE) ONLY. SERIOUS BODILY INJURY OR DEATH MAY RESULT IF A HAI SPEED HOLDER BROKE WHILE WORKING WITH SUPERCRITICAL FLUIDS.

DISCLAIMER

EXCEPT AS EXPRESSLY PROVIDED, ALL WARRANTIES, CONDITIONS, REPRESENTATIONS, INDEMNITIES AND GUARANTEES, WHETHER EXPRESS OR IMPLIED, ARISING BY LAW, CUSTOM, PRIOR ORAL OR WRITTEN STATEMENTS BY HIGGINS ANALYTICAL OR OTHERWISE (INCLUDING, BUT NOT LIMITED TO, ANY WARRANTY FOR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) ARE OVERRIDDEN, EXCLUDED AND DISCLAIMED.

EXCLUSIVE REMEDY

In the event a customer establishes a breach of the above warranty, customer's exclusive remedy against Higgins Analytical, and Higgins Analytical's sole liability, is, at its option, the replacement of the defective product or the refund of the purchase price therefor.

NO CONSEQUENTIAL DAMAGES
HIGGINS ANALYTICAL SHALL NOT BE RESPONSIBLE FOR ANY CONSEQUENTIAL, INDIRECT, SPECIAL, PUNITIVE OR INCIDENTAL DAMAGES, WHETHER FORESEEABLE OR UNFORESEEABLE, INCLUDING, BUT NOT LIMITED TO, CLAIMS FOR LOSS OF GOODWILL, PROFITS, USE OF MONEY OR USE OF PRODUCTS, STOPPAGE OF OTHER WORK OR IMPAIRMENT OF OTHER ASSETS, ARISING OUT OF BREACH OF EXPRESS OR IMPLIED WARRANTY, BREACH OF CONTRACT, MISREPRESENTATION, NEGLIGENCE, STRICT LIABILITY IN TORT OR OTHERWISE, EXCEPT ONLY IN THE CASE OF PERSONAL INJURY WHERE AND TO THE EXTENT THAT APPLICABLE LAW REQUIRES SUCH LIABILITY.

PATENTS AND TRADEMARKS

Kromasil is a trademark of Eka Nobel. NewGuard and MPLC are trademarks of Perkin-Elmer. PICOTAG is a trademark of Waters. HAI SIL, HAI CART, HAI PEEK, HAI FILTER, TARGA[®], CLIPEUS, PHALANX, HAI Speed, and HAIGUARD are trademarks of Higgins Analytical, Inc. Guard-Disc[®] is a trademark of Transgenomics, Inc. HAIGUARD and Guard-Disc[®] are covered by U.S. Patent 5,338,448 and other patents pending.

Higgins Analytical develops and manufactures a wide range of stainless steel and PEEK (polyetheretherketone) High Performance Liquid Chromatography (HPLC) columns and cartridges. These unique offerings address many needs of chromatographers desiring fast analysis times, efficient system and column protection, and chemical inertness.

Special application areas include LC/MS, combinatorial chemistry, column switching techniques, in-line filtration, and guard column requirements. Products are supplied and supported worldwide by a network of responsive distributors.

We are a quality-oriented chromatography company dedicated to providing performance products, expert support, and rapid delivery to research and quality control laboratories worldwide.



Higgins Analytical

*Expert Manufacturers of
HPLC Consumables*

Higgins Analytical, Inc.
1235 Pear Avenue, Suite 105
Mountain View, California, 94043 USA
Phone: 650-988-8930 or 800-988-6889
Fax: 650-988-8931 or 800-988-6881
E-mail hplc@higanalyt.com
www.higanalyt.com