

Direct Connect HPLC Columns

Applications:

- LC-MS
- Ballistic gradients
- High throughput screening

Unique Features:

- Tested and supplied with CofA
- PEEK/Titanium construction
- High performance
- Maximum short column efficiency
- Higgins Analytical Sorbents

Standard Sizes:


- 40 x 2.1mm
- 20 x 2.1mm
- 10 x 2.1mm

The *Sprite* family of narrowbore HPLC columns is a unique high performance product for LC-MS and other high-speed, small volume applications. All three sizes of 2.1mm I.D *Sprite* HPLC columns are individually tested. A *Certificate of Analysis* accompanies each column. As pictured to the right, the *Sprite's* certificate includes the original chromatogram and comprehensive testing and sorbent documentation.

The inert PEEK and titanium design has been proven to be structurally sound to over 8,000 PSI (540Bar).



Direct connect *Sprite* HPLC columns require no tools for leak-free finger-tight operation.



Typical *Sprite* Certificate of Analysis

Sprite Ordering Information



TARGA C18 5 μ m - *Sprite*
40 x 2.1mm
TR-0421-C185

20 x 2.1mm
TR-0221-C185

10 x 2.1mm
TR-0121-C185



CLIPEUS C18 5 μ m - *Sprite*
40 x 2.1mm
CR-0421-C185

20 x 2.1mm
CR-0221-C185

10 x 2.1mm
CR-0121-C185



HAISIL HL C18 5 μ m - *Sprite*
40 x 2.1mm
LR-0421-C185

20 x 2.1mm
LR-0221-C185

10 x 2.1mm
LR-0121-C185

For more information about Higgins Analytical Products, please ask us or your local distributor for our catalog. It can also be downloaded from our web site www.higanalyt.com



Higgins Analytical Products are sold and supported by distributors worldwide



Higgins Analytical, Inc.

Expert Manufacturer of HPLC Consumables

1235 Pear Ave, Suite 105, Mountain View, CA 94043

Phone: 650-988-8930 or 800-988-6889

Fax: 650-988-8931 or 800-988-6881

e-mail: hplc@higanalyt.com

Please visit our web site: www.higanalyt.com