



# CARBAMATE PESTICIDES

## HAISIL 100 C18



EPA 531.1

Column:  
HAISIL 100 C18 5µm  
250x4.6mm  
P/N HS-2546-M185




*n*-methylcarbamoyloximes and *n*-methylcarbamates in water  
Post column OPA derivatization and fluorescence detection

## CLYPEUS C18



EPA 632A

Column:  
CLYPEUS C18 5µm  
250x4.6mm  
P/N CS-2546-C185



Carbamates and Urea pesticides in industrial and municipal waste water, UV Detection